**IL-6**

**Lateral Flow Immunoassay for the semi-quantitative evaluation of human Interleukin-6 (IL-6)**  
English: Page 1-10

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

**Note:** Significant changes are indicated by dotted lines in the margin.  
A change history can be found at the end of the manual.

## Table of Contents

Explanation of Symbols .....	3
Warnings and Precautions .....	3
Materials Supplied, Storage and Stability .....	4
Materials Required (not included) .....	4
Intended Use .....	4
Specimen Collection and Preparation .....	4
Method and Test Principle .....	5
1) Test Performance with Milenia® POCScan Reader .....	6
2) Test Performance with Milenia® PicoScan .....	7
3) Test Performance with Evaluation Card (semi quantitative) .....	8
Interpretation of Results .....	9
Internal Quality Control .....	9
Assay Characteristics .....	10
Additional Products available .....	10

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## Warnings and Precautions

- Read **this instruction manual** carefully before performing the assay. Follow them step by step.
- Store reagents at 2 - 8° C.
- Do not use reagents after the expiry date has reached!
- Do not use Test Unit if foil pouch is damaged!
- Protect Test Units from humidity.
- Do not interchange components of packages of different batch codes.
- All samples, used Test Units, pipettes, tips, and other contaminated material should be handled carefully and as potentially infectious.
- The disposal of waste materials must be carried out according to current local regulations.
- For *in vitro* diagnostic use only. These reagents are to be used only by certified medical laboratory personnel authorized by the laboratory.
- The reagents provided in these kits contain materials of human and/or animal origin. Whenever human blood product is required for the preparation of these reagents, approved methods are used to test the blood for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.
- When participating in **RfB ring trials** (Reference Institute for Bioanalytics, Friesdorfer Str. 153, 53175 Bonn), the results must be entered under "**POCT- Point of Care Testing**".
- For professional users.

## Materials Supplied, Storage and Stability

Component	Cat.-No.	Content	Preparation	Storage	Shelf life
<b>Test Unit:</b> membrane coated with monoclonal (mouse) anti-hIL-6 antibody; monoclonal (mouse) anti-hIL-6 antibody conjugated to gold particles	MQSL6	20	ready to use	2 - 8 °C Use <b>immediately</b> after opening the foil pouch!	until expiry date
<b>Chase Buffer:</b> PBS-Buffer, pH 7.4	MQCB	1 vial à 2 mL	ready to use	2 - 8 °C	until expiry date
<b>Certificate of Analysis</b>	CoA-MQL6	1			
<b>Evaluation Card</b>	MQAL6	20			

Material Safety Data Sheets are available on request (see also [www.milenia-biotec.de](http://www.milenia-biotec.de)).

## Materials Required (not included)

- Pipet and tips for 100 µL.

## Intended Use

The Milenia® QuickLine IL-6 test is designed for the semi-quantitative evaluation of human Interleukin-6 (IL-6) in serum, plasma, cell culture supernatant, amniotic fluid or cerebrospinal fluid.

## Specimen Collection and Preparation

For determination of interleukin 6 (IL-6) serum, plasma (EDTA, heparinized, citrated), cell culture supernatant, amniotic fluid and cerebrospinal fluid can be used (**not whole blood**).  
Separate serum from blood cells within 4 hours. For storage longer than 8 hours at room temperature, samples should be aliquoted and stored frozen at -20 °C.  
The procedure calls for 100 µL sample volume per test.

## Method and Test Principle

The Milenia® QuickLine IL-6 test is a lateral flow immunoassay designed for the semi-quantitative evaluation of human Interleukin-6 (IL-6).

The sample is pipetted in the sample application port. IL-6 of patient's sample binds to a first monoclonal anti-IL-6 antibody conjugated to gold particles. The IL-6-loaded gold particles diffuse through the membrane and overflow the test line (T). There is a second monoclonal antibody specific for IL-6 coated on the membrane; so gold particles were bound specifically and become visible as a coloured line. Colour intensity is directly proportional to the concentration of IL-6 in the sample and intensifies during incubation time. The surplus of gold particles continues to diffuse over the test unit. The conjugate specific antibodies printed as a line on the membrane (control line, C) captures the rest of gold conjugate and a well visible line develops during the incubation time.

The IL-6 concentration correlates directly with colour's intensity of the test line (T).

Evaluation Options:

**1. Automatically:**

**a) Milenia® POCScan Reader**, semi-quantitative, page 6

The result will be generated by a stored lot specific standard curve.

For further information please read the Milenia® POCScan Reader Manual.

**b) Milenia® PicoScan**, semi-quantitative, page 7

The result will be generated by a stored standard curve. For

further information please read the Milenia® PicoScan Manual.



**2. Manually:**

Semi-quantitative, Evaluation Card (MQAL6), page 8

# 1) Test Performance with Milenia® POCScan Reader

**Important notes:**

- The test unit is sensitive to humidity. In case of damaged package do not use test unit!
- Before opening bring the test unit to room temperature (18 – 28 °C). Use the test unit immediately after opening. Do not store at room temperature!
- Prior performing the test read this instruction and the Milenia® POCScan Reader manual carefully! Follow the procedure step by step and notice the sequence!
- Start measurement immediately after pipetting the sample into the sample port! The incubation time of 20 minutes is controlled by the instrument.
- Before measurement, check the current lot of the used kit by outerbox label of the kit and cross check the kit lot with the *Certificate of Analysis* from this kit (look at manual: Milenia® POCScan Reader).


	Step	
1.		<p><b>Pipet 100 µL sample</b> in the round sample port of the test unit (see arrow in the figure).</p> <p><b>Attention:</b> Pipetting of <b>chase buffer</b> is <b>not</b> necessary!</p>
2.		<p>Close the drawer of the Milenia® POCScan Reader!</p> <p><b>Start measurement <u>immediately</u>.</b></p>
3.		<p>After <b>20 minutes</b> the <b>result</b> shows up in the display.</p> <p><b>The result must not be interpreted later!</b></p> <p><b>Only valid results are displayed.</b> If result is not valid, the test must be repeated. For details see chapter “Trouble Shooting” in the manual of the Milenia® POCScan Reader. The manual can be downloaded from the website <a href="http://www.milenia-biotec.de">www.milenia-biotec.de</a></p>

## 2) Test Performance with Milenia® PicoScan

### Important notes:


- The test is sensitive to humidity. In case of damaged package do not use test unit!
- Before opening bring the test unit and chase buffer to room temperature (18 – 28 °C). Use the test unit immediately after opening. Do not store at room temperature!
- Before performing the assay read this instruction and the PicoScan manual carefully! Follow the procedure step by step and notice the sequence!
- Start measurement immediately after adding the chase buffer to the patient sample in the sample port! The last incubation time of 5 minutes is controlled by the instrument.

1. Pipet **100 µL** sample into the sample application of the test unit (see figure).
2. After **15 minutes** incubation **add 2 drops chase buffer** directly with the bottle in the sample application of the test unit.
3. After further **5 minutes** of incubation, the test result is interpreted using the **PicoScan system**.

The result is automatically saved and can be directly read from the computer screen.


**The results should not be interpreted afterwards!**

### 3) Test Performance with Evaluation Card (semi quantitative)

**Important notes:**

- The test unit is sensitive to humidity. In case of damaged package do not use test unit!
- Before opening bring the test unit to room temperature (18 – 28 °C). Use the test unit immediately after opening. Do not store at room temperature!
- Before performing the assay read these instructions. Follow the procedure step by step and notice the sequence!
- The result must be interpreted exactly after 20 minutes.

	Step	
1.		<p><b>Pipet 100 µL sample</b> in the round sample port of the test unit (see arrow in the figure).</p>
2.		<p>After <b>15 minutes</b> incubation <b>add 2 drops chase buffer</b> directly with the bottle in the sample application of the test unit.</p>
3.	<p><b>Incubation</b></p>	<p>Incubate for <b>exactly 5 minutes</b></p>
4.		<p>After exactly <b>20 minutes</b>, the <b>result</b> must be interpreted with the provided Evaluation Card <b>The result must not be interpreted later!</b></p> <p>If result (see Internal Quality Control) is not valid, the test must be repeated.</p>

The IL-6 concentration of the sample correlates directly with colour's intensity of the test line (T). Evaluation of the Milenia® QuickLine IL-6 test is carried out by visual comparison of the test line colour intensity with the colour intensity of three different reference lines printed on the evaluation card provided with the kit.

The results are a semi-quantitative evaluation of the IL-6 concentration (pg/mL) present in the sample.

## Interpretation of Results

The IL-6 concentration of the sample correlates directly with colour's intensity of the test line (T).

### POCScan Reader/ PICOScan:

During measurement the software controls, whether the results are valid. If the test is not valid, no test result will be shown (look "Internal Quality Control").

Note: The devices cannot give an exact value. The values can only be regarded as approximated values.

### Evaluation Card:

Evaluation of the Milenia QuickLine IL-6 test is carried out by visual comparison of test line's colour intensity with the colour intensity of the three different reference lines printed on the evaluation card which is provided with the kit.

The results are a semi-quantitative evaluation of the IL-6 concentration (pg/mL) present in the sample.

The evaluation card may also be used for documentation. Beside the result, date of analysis, patient name, patient ID number, the operator and the lot number may be noticed on the card.


In General: Positive results should be verified concerning the entire clinical status of the patient. Additionally every decision for therapy should be taken individually.

## Internal Quality Control


**In any case, the control line (C) has to appear!**

It is used only as functional control and **cannot** be used for the interpretation of the test line. If the control line is not detectable after 20 minutes incubation time, the result is invalid! The analysis has to be repeated with a new test unit!

Valid results

Invalid results

## Assay Characteristics

Sample material:	Serum, plasma (EDTA, heparinized, citrated), culture supernatant, amniotic fluid, cerebrospinal fluid
Incubation time:	20 minutes at room temperature (18 - 28 °C)
Specificity:	Highly specific for human IL-6
Measurement range:	50 – 10,000 pg/mL
Sensitivity:	50 pg/mL
Reference values:	< 100 pg/mL in healthy blood donors
High-Dose-Hook-Effect:	Not detectable for samples ≤ 500,000 pg/mL IL-6
Interference Factors:	HAMA/Lipemia/Bilirubin/Haemolysis could interfere with this test. Elevated IL-6 level in such samples should be interpreted with caution
Calibration:	Calibrated with NIBSC WHO standard 89/548 (PicoScan, NIBSC IL-6)
Method Comparison:	60 plasma samples of ICU patients and 20 serum samples from healthy donors were analysed; Milenia® QuickLine IL-6 result were interpreted independently by two operators; 77.7 % (62/80) showed identical results in comparison to IMMULITE® IL-6 (DPC, Los Angeles, USA)

## Additional Products available

Product Name	Order No.	Content	Description
Milenia® QuickLine IL-6 Control	MQCO6 1	2 x 30 tests	The Controls are designed for internal quality control of Milenia® QuickLine IL-6 tests

For any questions and suggestions, please contact us:

**Telefon:** +49 641 – 94 88 83 – 0  
**Fax:** +49 641 – 94 88 83 – 80  
**Web:** [www.milenia-biotec.de](http://www.milenia-biotec.de)

**IL-6**

Lateral Flow Immunoassay for the semi-quantitative evaluation of human Interleukin-6 (IL-6)  
English: Page 1-10

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