Lactobacillus/ Pediococcus-Screen

Reagents for the detection of beer spoiling
Lactobacillus sp. and *Pediococcus* sp.
English: Page 3-11

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*** Attention:**

This kit contains only reagents used for amplification.
For the detection of amplicates the following additional kit is required:
Milenia® GenLine PCR-Universal module MGUP 1 (50 tests)

Table of Contents


Explanation of Symbols	3
Warnings and Precautions	3
Materials supplied, Storage and Stability	4
Materials Required (not included)	4
Intended Use	4
Method and Test Principle	4
Procedure of the Test – A short Overview	5
Preparation of Samples and PCR Reagents.....	6
Test Procedure.....	8
Interpretation of Results	9
Trouble shooting.....	10
Literature	10
Additional Products available.....	10

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Warnings and Precautions

- Read this instruction manual carefully before running the test. Follow the sequence of the work steps
- Pay attention to the expiration dates of the individual components of the test
- Complete thawing of all components prior use is required
- Do not exchange components with different Lot numbers
- Store all reagents at **-15 °C to -20 °C** in the original packaging
- **Immediately after use**, store all components at **-15 °C to -20 °C**
- The disposal of waste materials must be carried out according to current local regulations
- Wear laboratory coat and protective gloves
- In order to avoid contaminations, run the different steps of the test at different locations
- The place where the reaction mixes are prepared must be separated from the place where cultures are handled.
- Decontamination of the work area and the equipment used must be carried out regularly!

Materials supplied, Storage and Stability

REF: MGScLP 1						
Component	Part-Nr.	Quantity	Content	Preparation	Storage	Shelf life
PCR-Polymerase Master-Mix (2x)	MGPOL	2	300 µL	Ready to use	-15 °C to -20 °C	until expiry date
PCR- Primermix (2x)	MGPLP	2	280 µL	Ready to use	-20 °C	until expiry date
Positive Kit-Control	MGLHPC	1	250 µL	Ready to use	-20 °C	until expiry date
PCR-Water*	MGH2O	1	1 mL	Ready to use	-20 °C	until expiry date

*PCR-Water (MGH2O) can be used as diluent or negative control.

Material safety data sheets are available on request or from the company website (www.milenia-biotec.de)

Materials Required (not included)

Important:

In addition to the components supplied with the kit, the following reagents and items are required:

- **Milenia® GenLine PCR-Universal module (MGUP 1)**. Contains test strips and chase buffer.
- **Extraction Buffer (MGEXBA)**
- **Thermocycler Labcycler 48s**, Order code: MLCY, available at Milenia Biotec GmbH
- **Deep freezer: -20 °C**
- Pipets 2- 20 µL, 20- 200 µL, 100- 1,000 µL
- Sterile (DNA-free) pipet tips (with contamination protection for PCR)
- Sterile (DNA-free) reaction tubes 1.5 mL
- Sterile (DNA-free) PCR-tubes 0.2 mL e.g.: Supplier: Sorensen Bioscience Inc.
- 96-well microtiter plate
- (Desktop) centrifuge for reaction- and PCR-tubes
- Vortex mixer
- Magnetic stirrer
- Tube holder for reaction- and PCR-tubes, e.g.: isofreeze PCR tube rack

Intended Use

The Milenia® GenLine Lactobacillus/ Pediococcus-Screen is designed for the qualitative detection of beer spoiling organisms of the group *Lactobacillus* sp. and *Pediococcus* sp..

Method and Test Principle

Milenia® GenLine-Tests are rapid detection systems that combine PCR (Polymerase-Chain-Reaction) methodology with detection and interpretation of the resulting products by lateral flow tests.


In a first part of the test, labelled specific primers (short, single stranded DNA) are used to amplify specific DNA fragments. Thus the resulting PCR products carry the labels of the incorporated primers. In addition to the target gene, a control gene, which is also present in the PCR mixes, is amplified in order to make sure that the PCR process works properly.

In a second part of the test, the created PCR products are detected by a lateral flow Test Strip. The PCR products bind via their labels to antibodies attached to the solid phase of the strip and to colloidal gold particles which are coated with an antibody directed against the other label. Thus a “molecular sandwich” is formed and becomes visible as a line on the test strip.

Procedure of the Test – A short Overview

The entire procedure of the test can be divided into 3 different steps.

1. The PCR reagents and the samples are prepared.
2. After the addition of the sample to the PCR reagents, the PCR is started.
3. The resulting PCR products are detected by a simple lateral flow test strip.

Preparation of Samples and PCR Reagents

The preparation of samples and reagents includes the preparation of the PCR reaction mix and the samples to be analysed.

A) Preparation of PCR reaction mix:

Each individual tube of reaction mix has a total volume of 23 µL. In order to prepare the mix, the PCR-Polymerase Master Mix (2x) and the PCR-Primermix (2x) are mixed in a ratio as given in the table below. Finally the mix is filled in 23 µl aliquots into sterile PCR tubes.

The precise and contamination free preparation of the mix is the prerequisite of a proper function of the test. It is recommended to label all reaction tubes prior pipetting the samples. In case a positive (MGLHPC) and a negative (MGH2O) Control is included in the run, add 2 additional tubes to the number of samples introduced in the test run.

Number of PCR reaction tubes	PCR-Polymerase Master Mix (2x) (Volume to be pipetted)	PCR-Primermix (2x) (Volume to be pipetted)
1	13 µL	11 µL
5	65 µL	55 µL
10	130 µL	110 µL
20	260 µL	220 µL

Procedure:

- 1) According to the number of samples (plus a positive and negative control) PCR-Polymerase Master Mix (2x) is pipetted into a sterile 1.5 mL reaction tube (volume is given in the table).
- 2) PCR-Primermix (2x) is added to the tube (see table).
- 3) The mix of both components is mixed well (pipet, Vortex).
- 4) Take an appropriate number of sterile 0.2 mL PCR-reaction tubes.
- 5) Pipet 23 µL of mix into each PCR reaction tube.
- 6) Close the PCR-reaction tubes immediately and add the sample at a different location.

B) Preparation and addition of samples

No sample preparation is needed if the **analysis is done from liquid culture media**. The following media based on MRS can be used for the test:

Medium: NBB-B (Bouillon) – Döhler GmbH (Darmstadt)
 NBB-B-AM (Bouillon) – Döhler GmbH (Darmstadt)
 NBB-C (Concentrate) – Döhler GmbH (Darmstadt)

2µL from the enriched culture are transferred to the prepared PCR reaction mix.

Important Note: Very turbid enriched cultures contain large amounts of biological material. This can have a negative impact on the test or even may inhibit it completely. Therefore such samples must be pre- diluted 1:10 or 1:100 prior testing.

Analysis of single colonies from solid media (z.B.: BAT-A, Döhler GmbH):

Procedure:

- 1) Touch a colony on the agar plate with a pipet tip. Do not wipe across the colony.
- 2) Resuspend the cells in 100 µL sterile Water.
- 3) Add 2 µL of the cell suspension to the PCR reaction mix.

If too many cells are added to the mix, this may have a negative impact on the test. Therefore, a careful light touch in step 1) is recommended.

Test Procedure

The test procedure describes the start and termination of the PCR reaction as well as the subsequent lateral flow based detection of PCR products by a universal test strip (Milenia® GenLine Universal Module MGUP1).

C) PCR reaction

- 1) After addition of the sample to the PCR tube containing the reaction mix, close the tube and mix well.
- 2) Place all PCR tubes into a minicentrifuge and centrifuge the samples briefly.
- 3) All liquid should be free of air bubbles and located on the bottom of the PCR tubes.
- 4) Place all PCR tubes into a thermocycler.
- 5) Close the lid of the thermocycler.
- 6) Start the species-specific program within the folder „Bier“.
- 7) After completion of the program, stop the program and remove the PCR tubes from the thermocycler.
- 8) When the program is finished, the thermocycler holds a temperature of 8°C. PCR products must be stored in a refrigerator (app. 4°C) and should be analysed subsequently (max. 24 h later).

D) Detection of PCR products

For the detection of PCR products the **Milenia® GenLine PCR-Universal Module MGUP 1** is required. It is strongly recommended to use pipet tips with filter and separate pipets for the addition of the sample to be analyzed and the evaluation of PCR products to minimize the risk of DNA contamination. Regular decontamination of the pipets and a separation of the different working steps at different locations are recommended.

Procedure:

- 1) According to the number of tests, pipet 80 µL Chase buffer (MGCBA) each into clean wells of a 96 well microtiter plate.
- 2) Label each test strip (MGDS2A) needed for the analysis and take it aside.
- 3) Pipet 2 µL of each PCR product onto the sample application pad of a test strip.


- 4) Place the test strip with the sample application pad down into the Chase buffer of the wells of the microtiter plate and incubate for 5 min at room temperature. Read the results visually after exactly 5 min.
- 5) For the documentation of the results, strips can be fixed onto a data sheet (MGDBES). An additional documentation by taking an image of the strips with a camera is recommended.

Results must be read exactly after 5 minutes! After prolonged time unspecific, weak signals may occur on the test strip which are not associated with positive results.

Interpretation of Results

Tests are interpreted visually. Potential test results are shown in the figure below.

The control line [C] confirms the proper function of the test strip. The test line 1 [T1] is associated with the proper function of the PCR. The test line 2 [T2] is linked with the presence or absence of the target gene.

A) Positive Result (test strip 1, 2)

The prerequisite for a valid positive result is the appearance of a distinct test line [T1] (test strip 1 and 2).

In case very high levels of the target gene are introduced into the test, this leads to a diminished amplification of the amplification control. As a consequence no line may be visible at test line [T2]. Such results should be read positive (test strip 2).

In case of diminished signals (test strip 1* and 2*) it is recommended to rerun the samples after predilution.

B) Negative Result (test strip 3)

A valid negative result requires a distinct visible T2 line. This confirms that the PCR has worked properly. The test line 1 [T1] is not visible (test strip 3).

If the line of the PCR control [T2] is weakly visible, an inhibition of the PCR is very likely (test strip 3*). In this case it is recommended to rerun the samples after predilution.

C) No analysis possible (test strip 4, 5)

If neither the test line [T1] nor the PCR control [T2] is visible, no result can be reported (test strip 4). In case only the line at the location of the control line [C] is visible, either the PCR product has not been transferred successfully to the test strip or the amplification control has not performed properly.

If **no line** is visible on the test strip (test strip 5), please contact the customer service at Milenia Biotec GmbH

Trouble shooting

The negative control gives a positive result

DNA-Contaminations

Contaminations are a high risk when working with nucleic acids. It is very difficult to find out the reasons for contamination of materials of cellular origin or of purified DNA. Therefore the following preventive actions are strongly recommended:

- Check that the chosen PCR protocol of the thermocycler was correct.
- Decontaminate the workspace where the preparation of the reaction mix is done.
- Decontaminate the pipets.
- Use unopened, sterile plastic materials (Pipet tips/ 1.5 mL reaction vials/PCR vials)

After decontamination, a test run should be performed in order to check the effectiveness of the preventive actions.

It is recommended to run a positive and a negative control. The negative control **must** show a negative result.

If the negative control shows a positive result, this may be due to a contamination of the reagents. In this case new reagents (PCR-Water, PCR-Polymerase Master Mix (2x) and PCR-Primer Mix (2x) must be used. In the subsequent control run the negative control (PCR – Wasser, MGH2O) must be negative.

If the negative result shows a positive result again, please contact the customer service at Milenia Biotec GmbH.

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Koob, J: „A New Approach to Bacterial Identification in a Point-of-Care Format“; European Brewery Convention. Mundiconvenius, Porto, Portugal (2015)

Additional Products available

Product Name	Order No.	Content	Description
Milenia® <i>Lactobacillus brevis</i>	MGLBr Z*	24 tests	Reagents for the detection of <i>Lactobacillus brevis</i>
Milenia® <i>hop</i> resistance genes-Screen	MGSChOR 1*	48 tests	Reagents for the detection of <i>hop</i> resistance genes
Milenia® <i>Megasphaera/Pectinatus</i> -Screen	MGSsMP 1*	48 tests	Reagents for the detection of beer spoiling <i>Megasphaera sp.</i>

For any questions and suggestions, please contact us:

Telefon: 0 641 – 94 88 83 – 0
Fax: 0 641 – 94 88 83 – 80
E-mail: info@milenia-biotec.de
Web: www.milenia-biotec.de

Lactobacillus / Pediococcus-Screen

Reagents for the detection of beer spoiling *Lactobacillus* sp. and *Pediococcus* sp.

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Ist die negative Kontrolle erneut positiv, kontaktieren Sie den Kundenservice der Milenia Biotec GmbH.