## IL-6 Control

**Interleukin-6 quality control (bi-level) for  
Milenia® QuickLine IL-6 (MQL6 1)  
– Evaluation only with Milenia® POCScan Reader  
English: Page 1-2**

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


## Materials Supplied, Storage and Stability

Components	Cat. No.	Content	Preparation	Storage at	Shelf life
<b>QuickLine IL-6 Control 1:</b> Native, human IL-6 in serum matrix (concentration 1)	MQL6C1	3 vials à 1 mL	lyophilized; reconstitute with 1 mL aqua dest. <b>30 minutes</b> before use	2 - 8 °C; <b>Notice: After reconstitution, aliquot immediately the material and store at -20 °C!</b>	Until expiration date; after reconstitution 4 month at -20°C
<b>QuickLine IL-6 Control 2:</b> Native, human IL-6 in serum matrix (concentration 2)	MQL6C2	3 vials à 1 mL	lyophilized; reconstitute with 1 mL aqua dest. <b>30 minutes</b> before use	2 - 8 °C; <b>Notice: After reconstitution, aliquot immediately the material and store at -20 °C!</b>	Until expiration date; after reconstitution 4 month at -20°C

## Application

The controls are designed only for internal quality control of Milenia® QuickLine IL-6 tests (MQL6 1; test procedure see package insert). Test results can **only** be evaluated with Milenia® POCScan Reader.

The procedure calls for **100 µL** sample volume per test.

## Target Values

The values are test method-specific targets. Each laboratory should be establishing its own targets. Internal laboratory quality control must be adhering to current obligatory official regulations and directions.

Parameter	Cat. No.	Target (pg/ml)
QuickLine IL-6 Control 1	MQL6C1	<50
QuickLine IL-6 Control 2	MQL6C2	>300

## IL-6 Control

Interleukin-6 quality control (bi-level) for  
 Milenia® QuickLine IL-6 (MQL6 1)  
 – Evaluation only with Milenia® POCScan Reader  
 English: Page 1-2

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Milenia Biotec GmbH