HIT Positive Control

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Telefon: 0 641 – 94 88 83 – 0
 Fax: 0 641 – 94 88 83 - 80
www.milenia-biotec.de


HIT Positive Control

Qualitative *HIT Positive Control*, is exclusively designed for the application in Milenia® *QuickLine* HIT-Assay

REF:

MQCOHIT 1


Milenia Biotec GmbH
Versailler Str. 1
D-35394 Gießen
Germany

IVD

Fon.: +49- 6 41-94 88 83- 0
Fax: +49- 6 41-94 88 83- 80
E-mail: info@milenia-biotec.de
<http://www.milenia-biotec.de>

IVD

Note: Significant changes are indicated by dotted lines in the margin.
In the end of the IFU you will find a table with causes of changes.

MQCOHIT/B/2017/03/01

Intended Use

This reagent for internal quality control have to be used according to Milenia® *QuickLine* HIT-Assay's instructions of use to verify flow validity.

Kit Components, Storage and Shelf Life

Component	Cat.-No..	Content	Preparation	Storage	Shelf Life
QuickLine HIT Positive Control	MQHITC1	1 vial with lyophilisate, (40 tests)	Lyophilisate; reconstitute with 0,25 mL Aqua dest. not less than 30 min before usage	2°C to 8°C	lyophilised at 2°C to 8°C stable until the imprinted expiry date
Storage and components stability after reconstitution					
Control reconstituted in aqua dest.:				2 to 8 °C	1 month
Control is aliquoted immediately after reconstitution:				Freeze at -20 °C or below	3 months

Material safety data sheets are available on request.

Materials Required But Not Supplied

- Pipet and tips: 250 µL
- Aqua dest.
- Milenia® *QuickLine* HIT-Test (MQHIT 1 or MQHIT Z)
- Multi-axle rotating mixer

Sample volume

There is **5 µL Milenia® *QuickLine* HIT Positive Control** required per test unit (MQHIT 1 or MQHIT Z).

Preparation and Procedure

a) Please, follow these instructions for preparing control material for using the first time.

- Bring lyophilised Milenia® *QuickLine HIT Positive Control* to room temperature (18°C to 25°C) *before dissolving*.
- **Add 250 µL Aqua dest.**
- Close control vial.
- Mix the Milenia® *QuickLine HIT Positive Control* thoroughly with e.g. a Vortex Mixer. Afterwards put the control vial for at least 30 minutes, if possible, on a multi-axle rotating mixer.
- The Milenia® *QuickLine HIT Positive Control* is now ready-to-use.
- For long-term storage, it is recommended to aliquote (not required material) and freeze at, - 20 °C or below.

b) Milenia® *QuickLine HIT Positive Control*:

- **Add 5 µL** Milenia® *QuickLine HIT Positive Control* in the MQHIT test. Before test run please consult instructions for use for MQHIT 1 or MQHIT Z.
- Was the control material stored frozen or cooled, it should be brought to room temperature (at least 30 min on a multi-axle rotating mixer) before insert in the MQHIT Test.

Result

Milenia® *QuickLine HIT Positive Control* is evaluated positive, if the test line (T) shows a colored line, in a **darker and more intense shade of color, than the reference line on the evaluation card**.

Abb.: 1 - Positive

Internal quality control


Each laboratory should establish quality control procedures that are tailored to their individual requirements and are in compliant with relevant regulations or the terms of approval, including corrective actions which must be taken if the results are outside of target ranges.

In addition, the specifications in the instructions of Milenia® Quickline HIT tests (MQHIT 1 or MQHIT Z) is subject to the assessment of Milenia® Quickline HIT Positive Control:

Precautions

- Read these instructions for use and the instructions for use of Milenia® QuickLine HIT-Test (MQHIT 1/Z) carefully before performing the test.
- Follow them step-by-step
- Aliquoted control should be thawed completely before usage
- Bring reagents before using at room temperature
- Handle Milenia® QuickLine HIT Positive Control reagent carefully
- Wear protective clothing such as lab coats and disposable gloves
- Professional use
- Disposal must be carried out according to local regulations
- Animal materials which are used in the reagent are free from infection pathogens and are derived from animals kept in Germany.

Explanation of Symbols

Symbols	Explanation	Symbols	Explanation
	Expiry date		Consult Instructions for Use
	<i>In Vitro</i> Diagnostic Medical Device		Consult attended documents
	Batch code		Package size
REF	Catalogue number		Manufacturer
	Storage conditions		

Date/ Datum	Revision	Cause of Revision/ Änderungsgrund
28.02.2017	MQCOHIT 1_Z Rev B 2017-03-01	Registered Tradename / Handelsname changed / wird von Milenia® Quickline HIT Positive Kontrolle" in „Milenia® Quickline HIT Positive Control" geändert

For any questions and suggestions please contact us:

Telefon: 0 641 – 94 88 83 – 0
 Fax: 0 641 – 94 88 83 - 80
www.milenia-biotec.de