HIT Control

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

Note: Significant changes are indicated by dotted lines in the margin.
A change history can be found at the end of the manual.

Table of Contents

Explanation of Symbols	3
Warnings and Precautions	3
Materials supplied, Storage and Stability	4
Materials Required (not included)	4
Intended Use	4
Sample volume	4
Preparation and Procedure	5
Interpretation of Results	5
Internal quality control	6
Additional Products available	6

Explanation of Symbols

Symbols Symbole	Explanation (GB) Erklärung (DE)	Symbols Symbole	Explanation (GB) Erklärung (DE)
	Use by Verwendbar bis		Consult instructions for use Gebrauchsanweisungen beachten
	<i>In vitro</i> diagnostic medical device In-vitro Diagnostika		Consult attended documents Begleitdokumente beachten
	Batch code Chargenbezeichnung		Sufficient for <n> tests Ausreichend für <n> Prüfungen
	Catalogue number Katalognummer		Manufacturer Herstellr
	Storage conditions Lagerbedingungen		Control Kontrolle

Warnings and Precautions

- Read this instruction manual and the instruction manual of Milenia® *QuickLine* HIT-Test (MQHIT 1/Z) carefully before performing the assay
- Follow the sequence of the work steps
- Aliquoted control should be thawed completely before usage
- Allow all reagents reach room temperature before use
- Handle Milenia® *QuickLine HIT Control* reagents carefully
- Wear lab coat and protective gloves
- The disposal must be carried out according to local regulations
- Animal materials which are used in the reagent are free from infection pathogens and are derived from animals kept in Germany
- Human materials which are used in the reagent were tested and found negative for: HIV 1/2, HbsAg, HCV, HIV-1 (NAT) and RPR by FDA approved methods
- For professional users

Materials supplied, Storage and Stability

Component	Cat.-No.	Content	Preparation	Storage	Shelf Life
QuickLine HIT Positive Control	MQHITC1	1 vial with lyophilisate, (50 tests)	lyophilisate; reconstitute with 0,25 mL Aqua dest. not less than 30 min before usage	2°C to 8°C	lyophilised at 2°C to 8°C stable until the imprinted expiry date
QuickLine HIT Negative Control	MQHITC2	1 vial with lyophilisate, (50 tests)	lyophilisate; reconstitute with 0,25 mL Aqua dest. not less than 30 min before usage	2°C to 8°C	lyophilised at 2°C to 8°C stable until the imprinted expiry date

Storage and component stability after reconstitution	Storage	Shelf Life
Control reconstituted in aqua dest.	2°C to 8°C	1 month
Control is aliquoted immediately after reconstitution	-20°C or below	3 month

Safety data sheets are available on request (see also www.milenia-biotec.de)

Materials Required (not included)

- Pipet and tips: 250 µL
- Aqua dest.
- Milenia® *QuickLine* HIT-Test (MQHIT 1 or MQHIT Z)
- Multi-axle rotating mixer

Intended Use

The Milenia® *QuickLine HIT Control* is used for internal quality control of the corresponding Milenia® *QuickLine* HIT-Test. It is used to determine flow validity.

Sample volume

5 µL Milenia® *QuickLine HIT Positive/Negative Control* is required per test unit (MQHIT 1 or MQHIT Z).

Preparation and Procedure

a) Please follow these instructions to prepare the control material for first-time use

- Bring lyophilised Milenia® *QuickLine HIT Positive* and *Negative Control* to room temperature (18°C to 25°C) before dissolving
- **Add 250 µL aqua dest.** to each vial
- Close control vials
- Mix the Milenia® *QuickLine HIT Positive* and *Negative Control* thoroughly with e.g. a Vortex Mixer.
Afterwards put the control vials for at least 30 minutes on a multi-axle rotating mixer, if possible.
- The Milenia® *QuickLine HIT Control* is now ready-to-use.
- For long-term storage, it is recommended to aliquot (not required material) and freeze at - 20 °C or below.

b) Performance of Milenia® QuickLine HIT Control:

- **Add 5 µL** Milenia® *QuickLine HIT Positive* or *Negative Control* in the MQHIT test. Before run please consult instructions for use for MQHIT 1 or MQHIT Z.
- The control material should be brought to room temperature, if it was stored frozen or cooled (at least 30 min on a multi-axle rotating mixer) before using it in the MQHIT Test.

Interpretation of Results


Fig.1: Positive control

- Milenia® *QuickLine HIT Positive* Control is evaluated positive, if the test line (T) shows a colored line, in a **darker and more intense shade of color, than the reference line** on the evaluation card.


Fig.2: negative control

- Milenia® *QuickLine HIT Negative* Control is evaluated negative, if the test line (T) shows **no colored line**.

Internal quality control

Each laboratory should establish quality control procedures that are tailored to their individual requirements and are in compliance with relevant regulations or the terms of approval, including corrective actions which must be taken if the results are outside of target ranges.

In addition, the specifications in the instructions of Milenia® QuickLine HIT tests (MQHIT 1 or MQHIT Z) is subject to the evaluation of Milenia® Quickline HIT Control

Additional Products available

Product Name	Order No.	Content	Description
Milenia® QuickLine HIT	MQHIT 1	20 tests	Qualitative detection of IgG antibodies against PF4/polyanion-complexes in human citrated plasma or serum
	MQHIT Z	5 tests	

For any questions and suggestions, please contact us:

Telefon: +49 641 – 94 88 83 – 0
Fax: +49 641 – 94 88 83 – 80
Web: www.milenia-biotec.de

HIT Control

Qualitative *HIT Control*, exclusively designed for the application in Milenia® *QuickLine* HIT-Assay

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