

Detection of PCR-Products

Comparison of the Agarose Gel and Milenia® HybriDetect

	Agarose Gel	Milenia® HybriDetect
Time	Preparation Preparation of Running Buffer 5 min Weigh and Boil Agarose 5 min Allow Agarose to solidify app. 30 min Overload with Gel Buffer 1 min Load Gel app. 5 min	Label Strip Pipet Running Buffer app. 5 min Transfer Sample to Strip
	Incubation Time Run Gel app. 1 h	Run Strip 5 min
	Visualization Visualization app. 10 min	Visualization app. 2 min
Total Time	app. 1 to 2 h	app. 15 min
Additional Items Required	Instruments Pipets, Scale, Microwave, Power Supply, Electrophoresis Chamber, Gel-Visualization Device, UV-Table	Pipets, Pen, Camera / Smart Phone
	Reagents Running Buffer (TBE or TAE) Loading Buffer Ethidium bromide	Strips (MGHD1/MGHD2) Running Buffer (MGCB/MGCB2)
	Sources of Danger Electric Power UV-Light Ethidium bromide	none
	Total	many
Price	Instruments High Investment Costs	Low Investment Costs
	Primers Unlabeled Primers → Low Costs	Labeled Primers → Higher Costs (app. + 30 € per Primer)
	Reagents Agarose, Running Buffer, Loading Buffer, Molecular Weight Standard and Ethidium bromide or expensive DNA-Stainers (Ethidium bromide (Waste): Costs for Disposal) → app. 0,90 €/small Gel (10 Slots)	→ app. 2,50 € per Strip
	Personnel Costs Educated Personnel necessary High Investment in Time	No Educated Personnel necessary Low Investment in Time
	Total Costs	low (if Sample Throughput is high)
Use	Experience / Education necessary	not necessary
	Sensitivity high	very high (app. 100 - 1000 x higher)
	Selectivity good	very good (Distinguishes between PCR-Products of the same size!)
	Application Many fold	restricted
	Distribution (Market) high	low

