milenia biotec

Milenia QuickLine HAMA

English

Lateral Flow Immunoassay to detect human anti-mouse antibodies (HAMA)

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
Σ 5

MQHM Z / B / 2023-05-19

Note: Significant changes are indicated by dotted lines in the margin.
A change history can be found at the end of the manual.

Milenia Biotec GmbH

Versailler Str. 1
D-35394 Gießen
Germany

Tel.: +49-641-94 8883-0

Fax: +49-641-94 8883-80

E-mail: info@milenia-biotec.de

Web: www.milenia-biotec.com

IVD


CE

Table of Content

English

Explanation of Symbols	3
Warnings and Precautions	3
Materials Supplied, Storage and Stability	4
Materials Required but not Supplied	4
Intended Use	4
Specimen Collection and Preparation	5
Volume Required	5
Summary and Explanation	5
Test Principle	5
Test Procedure	6
Interpretation of Results	7
Function Control	8
Assay Characteristics	8
Limitation/ Interference Factors	8
Contact	9
Change History	9

Explanation of Symbols

Symbol	Explanation	Symbol	Explanation
	Manufacturer		Do not reuse
	Use by		Consult Instructions for Use
	Batch Code		Consult attended documents
	Catalogue number		<i>In-vitro</i> diagnostic medical device
	Do not use if package is damaged		Contains sufficient for <n> tests
	Storage conditions		CE Marking

Warnings and Precautions

- Read **this instructions for use** carefully before performing the test and follow step by step.
- Store reagents at 2 - 8° C.
- Do not use Test Unit if foil pouch is damaged.
- Protect Test Units from humidity.
- Do not interchange components of packages of different batch codes.
- All samples, used Test Units, pipettes, tips, and other contaminated material should be handled carefully and as potentially infectious.
- The disposal of waste materials must be carried out according to current local regulations.
- For *in vitro* diagnostic use only.
- These reagents are to be used only by certified medical laboratory personnel authorized by the laboratory.
- The reagents provided in these kits contain materials of human and/or animal origin. Whenever human blood product is required for the preparation of these reagents, approved methods are used to test the blood for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.
- For professional users

Materials Supplied, Storage and Stability

Materials Supplied

REF	Contents	Number of Tests
MQHM Z	5 x HAMA Test Units (MQSHM) 1 x 0,6 mL Chase Buffer (MQBHM) 5 x Evaluation Cards (MQAHMZ)	5

Storage and Stability

Components	Description	Preparation	Storage	Shelf Life
HAMA Test Units	Membrane coated with mouse IgG; Mouse IgG conjugated to gold nanoparticles	Ready to use	2-8°C Use immediately after opening the foil pouch!	Until expiry date In the closed foil pouch!
Chase Buffer	pH 7.2	Ready to use	2-8°C	Until expiry date
Evaluation Card	Result Interpretation only by this card via colored line	Ready to use	Protected from light, in the associated test kit	

Material Safety Data Sheets are available on request, alternatively on the relevant product page on our website (please see www.milenia-biotec.de)

Materials Required but not Supplied

- Pipet and tips: 50 µL

Intended Use

- For *in vitro* diagnostic use.
- The HAMA kit is designed for the detection of human anti-mouse antibodies (IgG).

Specimen Collection and Preparation

Serum

Use tubes without any anticoagulant!

Please Note!

- Avoid haemolysis
- Separate serum from blood cells within two hours.
- If there is no centrifuge available, place serum sample tubes after blood collection (as described in sample tubes' IFU) vertically in a sample rack until the natural sedimentation has finished and a stable blood clot has been formed.

Volume Required

- 50 µl serum

Summary and Explanation

Monoclonal mouse-antibodies are used in vivo e.g. to prevent rejection of the graft or for immunotherapy in tumour patients, respectively. During the course of such therapy patient's can form antibodies against the mouse-antibodies (HAMAs). The formation of HAMAs can have an impact on the effectiveness of the therapy. In addition, HAMAs lead to troublesome cross-reactions in immunoassays, which can cause incorrect results.

Test Principle

The HAMA is a lateral flow immunoassay designed for qualitative determination of human anti-mouse antibodies in human serum.

The sample will be pipetted in the sample application port of the test unit and the chase buffer is added immediately. The buffer forces the sample to migrate through the membrane of the test unit.

Human anti-mouse antibodies from patient's sample bind first to mouse antibodies which are conjugated to gold nanoparticles.

In a second step the human anti-mouse antibodies- from patient's sample- will be caught by mouse IgG antibodies, coated on the T-line of the nitrocellulose membrane.

In the presence of the HAMA-Conjugat-complexes a band becomes visible at the location of the test line (T). This test line- if it is positive- has to be the same or a darker shade of color than the given one on the provided evaluation card.


The surplus of gold particles continues to migrate through the membrane and is captured at the control line (C) by specific antibodies. A well visible line appears after incubation time.

Color's intensity of the test line is directly proportional to the HAMA concentration in the sample and increases during incubation time.

Test Procedure


General Remarks:

- Before opening bring test unit and sample buffer to room temperature (18 – 28 °C)!
- Use test unit immediately after opening.
- Do not store test unit at room temperature!

<p>1. Sample</p>	<p>Add 50 µL serum to the round sample port of the test unit</p>	
<p>2. Buffer</p>	<p>Add 1 drop chase buffer</p> <p>Note: Do not contaminate the vial with sample material! Avoid air bubbles!</p>	
<p>3. Incubation</p>	<p>10 minutes at room temperature (18 - 28 °C)</p>	
<p>4. Evaluation</p>	<p>The result can only be interpreted visually by the provided Evaluation Card: The test is interpreted positive if the test line shows the same or a darker shade of color after incubation time than the given line in the evaluation card (for more information see Interpretation of Results and Functional Control). The result must be interpreted immediately after incubation time.</p>	

Interpretation of Results

Hold the colored line of the Evaluation Card (as shown below) close beside the test line of your test unit (See figures 1, 2).

Negative Results

- **Test Line (T)** shows a **lighter shade of colour** than the given line on the Evaluation Card
- **No Test Line (T)** is visible

Positive Results

- **Test Line (T)** shows the **same (see Fig. 1) or a darker (see Fig. 2) shade of colour** than the given line on the Evaluation Card

The **Control Line-C** must be visible for a valid result. (see **Function Control** of the test unit)

Please note:

Results should be verified concerning the entire clinical status of the patient.
Additionally every decision for therapy should be taken individually.


Function Control

In any case, the Control line (C) has to appear!


This line is only used as functional control and is no result.

If the control line is not detectable after 10 minutes incubation time, the result is invalid and the analysis has to be repeated with a new test unit!

Valid results

Invalid results

Assay Characteristics

Sample material: 50 μ L serum

Incubation time: 10 minutes at room temperature (18 – 28 °C)

Diagnostic sensitivity: 97%

Diagnostic specificity: 91%

Limitation/ Interference Factors

Do not use plasma samples!

Contact

For additional information and support, please contact us:

**Milenia Biotec GmbH**

Versailler Str. 1

D-35394 Gießen

Germany

Tel.: +49-641-94 8883-0

Fax: +49-641-94 8883-80

E-mail: info@milenia-biotec.deWeb: www.milenia-biotec.com

Change History

Date	Revision	Cause of Revision
23/04/2013	A	Creation of the HAMA instructions for use
19/05/2023	B	Formal revision

milenia biotec

Milenia QuickLine HAMA

Deutsch

Lateral Flow Immunoassay to detect human anti-mouse antibodies (HAMA)
English: Page 1 - 9

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Milenia Biotec GmbH · Versailler Str. 1 · 35394 Gießen, Germany