milenia biotec


IL-6

**Lateral Flow Immunoassay for the semi-quantitative
evaluation of human Interleukin-6 (IL-6)**
English: Page 1-7

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E-Mail: info@milenia-biotec.de

Web: www.milenia-biotec.com

Note: Significant changes are indicated by dotted lines in the margin.
A change history can be found at the end of the manual.

IFU / REF ML6 1 / REV ZC / 2024-10-15

Table of Contents

Explanation of Symbols	3
Warnings and Precautions	3
Materials Supplied, Storage and Stability	4
Materials Required (not included)	4
Intended Use	4
Specimen Collection and Preparation	4
Method and Test Principle	4
Test Performance with Evaluation Card	5
Interpretation of Results	6
Internal Quality Control	6
Assay Characteristics	7
References	7

Explanation of Symbols

Symbol	Explanation	Symbols	Explanation
	Use-by date		Contains sufficient for <n> tests
	<i>In vitro</i> diagnostic medical device		Manufacturer
	Batch code		Do not use if package is damaged and consult instructions for use
	Catalogue number		Consult instructions for use
	Temperature limit		Contains biological material of animal origin
	Do not re-use		Unique Device Identifier

Warnings and Precautions

- Read **this instruction manual** carefully before performing the assay. Follow them step by step.
- Store reagents at 2 - 8° C.
- Do not use reagents after the expiry date has reached!
- Do not use Test Unit if foil pouch is damaged!
- Protect Test Units from humidity.
- Do not interchange components of packages of different batch codes.
- All samples, used Test Units, pipettes, tips, and other contaminated material should be handled carefully and as potentially infectious.
- The disposal of waste materials must be carried out according to current local regulations.
- For *in vitro* diagnostic use only. These reagents are to be used only by certified medical laboratory personnel authorized by the laboratory.
- The reagents provided in these kits contain materials of human and/or animal origin. Whenever human blood product is required for the preparation of these reagents, approved methods are used to test the blood for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.
- When participating in **RfB ring trials** (Reference Institute for Bioanalytics, Friesdorfer Str. 153, 53175 Bonn), the results must be entered under "**Interleukin-6-Streifentest**".
- For professional users.

Materials Supplied, Storage and Stability

Component	Cat.-No.	Content	Preparation	Storage	Shelf life
Test Unit: membrane coated with monoclonal (mouse) anti-hIL-6 antibody; monoclonal (mouse) anti-hIL-6 antibody conjugated to gold particles	MQSL6	20	ready to use	2 - 8 °C	until expiry date. Use immediately after opening the foil pouch!
Chase Buffer: PBS-Buffer, pH 7.4	MQCB	1 vial à 2 mL	ready to use	2 - 8 °C	until expiry date
Evaluation Card	MQAL6	20			

Material Safety Data Sheets are available on request (see also www.milenia-biotec.com).

Materials Required (not included)

- Pipet and tips for 100 µL.

Intended Use

The Milenia QuickLine IL-6 test is an *in vitro* diagnostic lateral flow immunoassay for the semiquantitative detection of human Interleukin-6 in serum, plasma, or amniotic fluid.

Specimen Collection and Preparation

For determination of Interleukin-6 (IL-6) serum (2), plasma (EDTA, heparinized, citrated) and amniotic fluid (1) can be used (**not whole blood**).

Separate serum from blood cells within 4 hours. For storage longer than 8 hours, samples should be aliquoted and stored frozen at -20 °C.

The procedure requires for 100 µL sample volume per test.

Method and Test Principle

The Milenia QuickLine IL-6 test is a lateral flow immunoassay designed for the semi-quantitative evaluation of human Interleukin-6 (IL-6).

The sample is pipetted in the sample application port of the test unit. IL-6 of patient's sample binds to a first monoclonal anti-IL-6 antibody conjugated to gold particles. The IL-6-loaded gold particles diffuse through the membrane and overflow the test line (T). There is a second monoclonal antibody specific for IL-6 coated on the membrane; so gold particles were bound specifically and become visible as a coloured line. Colour intensity is directly proportional to the concentration of IL-6 in the sample and intensifies during incubation time. The surplus of gold particles continues to diffuse over the test unit. The conjugate specific antibodies sprayed as a line on the membrane (control line, C) captures the rest of gold conjugate and a well visible line develops during the incubation time.


The IL-6 concentration correlates directly with colour's intensity of the test line (T).

The semi-quantitative evaluation is carried out manually using the provided Evaluation Card (MQAL6) (see page 5).

Test Performance with Evaluation Card

Important notes:

- The test unit is sensitive to humidity. In case of damaged package do not use test unit!
- Before opening bring the test unit to room temperature (18 – 28 °C). Use the test unit immediately after opening. Do not store at room temperature!
- Before performing the assay read these instructions. Follow the procedure step by step and notice the sequence!
- The result must be interpreted exactly after 20 minutes.

	Step	
1.		Pipet 100 µL sample in the round sample application port of the test unit (as shown in the picture).
2.		After 15 minutes incubation add 2 drops chase buffer directly with the bottle in the round sample application port of the test unit.
3.	Incubation	Incubate for exactly 5 minutes
4.		After exactly 20 minutes , the result must be interpreted with the provided Evaluation Card The result must not be interpreted later! If result (see Internal Quality Control) is not valid, the test must be repeated.

The IL-6 concentration of the sample correlates directly with colour's intensity of the test line (T). Evaluation of the Milenia QuickLine IL-6 test is carried out by visual comparison of the test line colour intensity with the colour intensity of three different reference lines printed on the Evaluation Card provided with the kit.

The results are a semi-quantitative evaluation of the IL-6 concentration (pg/mL) present in the sample.

Interpretation of Results

The IL-6 concentration of the sample correlates directly with colour's intensity of the test line (T).

Evaluation of the Milenia QuickLine IL-6 test is carried out by visual comparison of test line's colour intensity with the colour intensity of the three different reference lines printed on the Evaluation Card which is provided with the kit.

The results are a semi-quantitative evaluation of the IL-6 concentration (pg/mL) present in the sample.

The Evaluation Card may also be used for documentation. Beside the result, date of analysis, patient name, patient ID number, the operator and the lot number may be noted on the card.


Note: Positive results should be verified concerning the entire clinical status of the patient. In addition, every decision for therapy should be taken individually.

Internal Quality Control


In any case, the control line (C) must appear!

The C-line is used only as functional control and **cannot** be used for the interpretation of the test line. If the control line is not detectable after 20 minutes incubation time, the result is invalid! The analysis has to be repeated with a new test unit!

Valid results

Invalid results

Assay Characteristics

Sample material:	Serum, plasma (EDTA, heparinized, citrated), amniotic fluid
Incubation time:	20 minutes at room temperature (18 - 28 °C)
Specificity:	Amniotic fluid: 96% (1) Serum: 82% (2)
Measurement range:	100- 10,000 pg/mL (Evaluation Card)
Sensitivity:	Amniotic fluid: 97% (1) Serum: 100% (2)
Analytical Sensitivity:	100 pg/mL (Evaluation Card)
Reference values:	< 100 pg/mL in healthy blood donors
High-Dose-Hook-Effect:	Not detectable for samples $\leq$ 500,000 pg/mL IL-6
Interference Factors:	HAMA/Lipemia/Bilirubin/Haemolysis could interfere with this test. Elevated IL-6 level in such samples should be interpreted with caution
Calibration:	Calibrated with NIBSC WHO standard 89/548
Method Comparison:	60 plasma samples of ICU patients and 20 serum samples from healthy donors were analysed; Milenia QuickLine IL-6 result were interpreted independently by two operators; 77.5 % (62/80) showed identical results in comparison to IMMULITE® IL-6 (DPC, Los Angeles, USA)

References

1. Musilova, I. et al. Intraamniotic inflammation in women with preterm prelabor rupture of membranes. PLoS One 10, (2015).
2. Schlosser, H.-G., Brock, M., Woiciechowsky, C. & Volk, H. Interleukin 6 Level Determination in Severe Head Injury Using Bedside Densitometry. Crit. Care Med. 33, A47 (2005).

For any questions and suggestions, please contact us:

Phone : +49 641 – 94 88 83 – 0
Fax: +49 641 – 94 88 83 – 80
Web: www.milenia-biotec.com

IL-6

Lateral Flow Immunoassay for the semi-quantitative
evaluation of human Interleukin-6 (IL-6)
English: Page 1-7

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Milenia Biotec GmbH · Versailler Str. 1 · 35394 Gießen, Germany