Milenia GenLine

PCR Universal-Modul

Deutsch

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DE: Seite 1-11

Lateral flow dipsticks for the detection of genomic
amplificates of beer spoiling germs.

EN: Page 11-22

REF MGUP 1*

 50

IFU / REF MGUP 1 / REV A / 2025-07-14

*** Attention:**

This kit (Universal Module) only contains reagents for the detection of the
amplificate.

An additional, specific test kit (PCR Module) is required for the previous PCR
reaction.

Note: Significant changes are indicated by dotted lines in the margin.
A change history can be found at the end of the instructions for use.

Milenia Biotec GmbH

Versailler Str. 1

35394 Gießen

Germany

Tel.: +49 641 94 8883-0

Fax: +49 641 94 8883-80

E-mail: info@milenia-biotec.de

Web: www.milenia-biotec.com

Table of Contents

English	12
Explanation of Symbols.....	14
Warnings and Precautions	14
Materials Supplied, Storage and Stability	15
Safety Data Sheets (SDS)	16
Materials Required but not Supplied	16
Intended Use.....	16
Method and Test Principle	16
Procedure of the Test – A short Overview.....	17
Step 1: Preparation of PCR Reagents and Samples /	17
Step 2: PCR Reaction.....	17
Step 3: Detection of PCR Products - Evaluation by a simple lateral flow test strip	18
Interpretation of Results	19
a) Positive Result	19
b) Negative Result	19
c) No Analysis Possible	20
Troubleshooting.....	20
References	20
Additional Products Available	21
Contact.....	21
Change History	22

Explanation of Symbols

Symbol	Explanation	Symbol	Explanation
	<i>Manufacturer</i>		Contains sufficient for <n> tests
	Use-by date		Consult <i>instructions for use</i> or consult electronic <i>instructions for use</i>
	<i>Batch code</i>		Unique device identifier
	<i>Catalogue number</i>		Temperature limit
	Contains biological material of animal origin		

Warnings and Precautions

- Store all reagents at 2 - 8° C in their original containers.
- Before use, bring all reagents to room temperature (18 – 28° C).
- The expiration date of all components must be observed. Do not use components past the expiration date!
- Protect Dipsticks from humidity. Containers must always be closed.
- Touch and label only the foil- covered areas of the dipsticks.
- Do not interchange components from different lots.
- The disposal of waste materials must be carried out according to current local regulations.
- The assay buffer contains an anti-microbial reagent; therefore, avoid contact with skin and/or mucous membranes.
- Animal materials which are used in the reagent are free from infection pathogens and are derived from animals kept in Germany
- For professional users

Materials Supplied, Storage and Stability

Materials Supplied

REF	Content	Part.-No.	No. of Tests
MGUP 1	50	Dipsticks for Two Test Bands	50
	10 mL	Chase Buffer	

Storage and Stability

Kit-Component	Description	Preparation	Storage	Shelf Life
HybriDetect Dipsticks (MGDS2A)	Membrane coated with biotin-ligand, polyclonal anti-digoxigenin antibody (goat) and polyclonal anti-goat antibody (rabbit) Gold coated with polyclonal anti-FAM/FITC antibody (goat)	ready-to-use, bring to room temperature	2-8 °C Always keep the container closed to protect it from moisture!	until expiry date
Chase Buffer (MGCBA)	Tris-buffered Saline	ready-to-use, bring to room temperature	2-8 °C	until expiry date

Safety Data Sheets (SDS)

Safety Data Sheets (SDS) are available on request and alternatively under the respective product on our website (www.milenia-biotec.com).

[SDS REF MGUP 1 - Universal Module](#)

Materials Required but not Supplied

- Pipettes
- Pipette tips (containing protective filters for PCR)
- Reaction tubes or 96-well microtiter plate

Intended Use

The PCR Universal Module consists of lateral flow dipsticks and is used for the detection of Genomic Amplificates (beer germs).

Method and Test Principle

Milenia GenLine-Tests are rapid detection systems that combine PCR (Polymerase-Chain-Reaction) methodology with detection and interpretation of the resulting products by lateral flow tests.

This kit includes all components needed for the **evaluation of the PCR tests** and is therefore the **second part of the entire test**.

In a first part of the test, labelled specific primers (short, single stranded DNA) are used to amplify specific DNA fragments. Thus, the resulting PCR products carry the labels of the incorporated primers. In addition to the target gene, a control gene, which is also present in the PCR mixes, is amplified to make sure that the PCR process works properly.

In a second part of the test, the created PCR products are detected by a lateral flow test strip. The PCR products bind via their labels to antibodies attached to the solid phase of the strip and to colloidal gold particles which are coated with an antibody directed against the other label. Thus a "molecular sandwich" is formed and becomes visible as a line on the test strip.

Procedure of the Test – A short Overview

The entire procedure of the test can be divided into 3 different steps.

1. Preparation of PCR reagents and samples (**PCR Module**)
2. PCR reaction (**PCR Module**)
3. Evaluation by a simple lateral flow test strip (**PCR Universal Module**)


Figure 2: Presentation of how the procedure is carried out

Step 1: Preparation of PCR Reagents and Samples /

Step 2: PCR Reaction

The preparation of samples and reagents includes the preparation of the PCR reaction mix and the samples to be analysed. For this, a Milenia GenLine-Test (PCR Module) suitable for the specific application is required. The tests include the PCR-Polymerase Master-Mix, PCR-Primermix, Positive Kit-Control and PCR-Water.

For more information regarding the preparation of samples and reagents as well as the PCR procedure, please refer to the instructions for use of the required Milenia GenLine-Test.

Available PCR Modules:

[Milenia GenLine Hop Resistance Gene Screen \(REF: MGScHOR 1\)](#)

[Milenia GenLine *Lactobacillus/Pediococcus* Screen \(REF: MGScLP 1\)](#)

[Milenia GenLine *Megasphaera/Pectinatus* Screen \(REF: MGScMP 1\)](#)

Further information on available products can be found in the chapter [Additional Products Available](#).

Step 3: Detection of PCR Products - Evaluation by a simple lateral flow test strip

Procedure:

1. According to the number of tests, pipet 80 μ L chase buffer (MGCBA) each into clean wells of a 96-well microtiter plate.
2. Label each test strip (MGDS2A) needed for the analysis and take it aside.
3. Pipet 2 μ L of each PCR product onto the sample application pad of a test strip.

4. Place the test strip with the sample application pad down into the chase buffer of the wells of the microtiter plate and incubate for 5 min at room temperature. Read the results visually after exactly 5 min.
5. A documentation by taking an image of the strips with a camera is recommended.

Note 1:

Results must be read exactly after 5 minutes! After prolonged time unspecific, weak signals may occur on the test strip which are not associated with positive results.

Note 2:

For the detection of PCR products it is strongly recommended to use pipet tips with filter and separate pipets for the addition of the sample to be analyzed and the evaluation of PCR products to minimize the risk of DNA contamination. Regular decontamination of the pipets and the separation of the different working steps at different locations are recommended.

Interpretation of Results

Tests are interpreted visually. Potential test results are shown in the figure below.


Figure 2: Potential test results for visual interpretation

The control line [C] confirms the proper function of the test strip. The test line 1 [T1] is associated with the presence or absence of the target gene. The test line 2 [T2] is linked with the proper function of the PCR.

a) Positive Result

The prerequisite for a valid positive result is the appearance of a distinct test line [T1] (test strip **1 and 2**). If high levels of the target gene are introduced into the test, the amplification control can vanish. Nevertheless, there will be a strong signal on test line T1. Therefore, no line may be visible at the test line [T2]. Such results should be read positive (test strip **2**).

b) Negative Result

A valid negative result requires a distinct visible test line [T2]. This confirms that the PCR has worked properly. The test line [T1] is not visible (test strip **3**).

Note:

If the line of the PCR control [T2] is weakly visible, an inhibition of the PCR is very likely (test strip **3***). In this case it is recommended to rerun the samples after predilution.

c) No Analysis Possible

If neither the test line [T1] nor the PCR control [T2] is visible, no result can be reported (test strip **4**). In case only the line at the location of the control line [C] is visible, either the PCR product has not been transferred successfully to the strip or the amplification control has not performed properly.

If **no line** is visible on the test strip (test strip **5**), please contact the customer service at Milenia Biotec GmbH.

Troubleshooting

The negative control gives a positive result.

Contaminations are a high risk when working with nucleic acids. It is very difficult to find out the reasons for contamination of materials of cellular origin or of purified DNA. Therefore, the following preventive actions are strongly recommended:

- Check that the chosen PCR protocol of the thermocycler was correct.
- Decontaminate the workspace where the preparation of the reaction mix is done.
- Decontaminate the pipets.
- Use unopened, sterile plastic materials (pipet tips/ 1.5 mL reaction vials/PCR vials).

After decontamination, a test run should be performed in order to check the effectiveness of the preventive actions. It is recommended to run a positive and a negative control. The negative control **must** show a negative result.

If the negative control shows a positive result, this may be due to contamination of the reagents. In this case new reagents (PCR-Water, PCR-Polymerase master mix (2x) and PCR-Primer mix (2x) must be used. In the subsequent control run the negative control (PCR-Water, MGH2O) must be negative.

If the negative result shows a positive result again, please contact customer service at Milenia Biotec GmbH.

References

1. Koob, J et al., „Molekularbiologischer Nachweis bierschädigender Organismen, Der Weihenstephaner, October 2015
2. Koob, J: „A New Approach to Bacterial Identification in a Point-of-Care Format“; European Brewery Convention. Mundiconvenius, Porto, Portugal (2015)

Additional Products Available

REF	Produktname	Testanzahl	Beschreibung
MGSMP 1	<u>Megasphaera/ Pectinatus-Screen</u>	48 tests	Reagents for the detection of beer spoiling <i>Megasphaera sp.</i> and <i>Pectinatus sp.</i>
MGSCHOR 1	<u>Hop Resistance Gene-Screen</u>	48 tests	Reagents for the detection of hop resistance genes.
MGSCLP 1	<u>Lactobacillus/ Pediococcus-Screen</u>	48 tests	Reagents for the detection of beer spoiling <i>Lactobacillus sp.</i> and <i>Pediococcus sp.</i>
MGES 1	<u>Extraction System</u>	1000 reactions	Universal kit for the isolation of analytes from swabs for further processing.
MGSDB	<u>Extraction Buffer</u>	100 reactions	Buffer for purification of swab samples with MGES 1.

Contact

For additional information and support, please contact us:

Milenia Biotec GmbH

Versailler Str. 1

35394 Gießen

Germany

Tel.: +49 641 94 8883-0

Fax: +49 641 94 8883-80

E-mail: info@milenia-biotec.de

Web: www.milenia-biotec.com

Change History

Date	Revision	Cause of Revision
2025-07-14	IFU/REF MGUP 1/REV A / 2025-07-14	IFU created