Milenia GenLine **Extraktionssystem**

Deutsch

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DE: Seite 1-7

Extraction system for detaching
sample material from swabs

EN: Page 8-14

REF MGES 1

 1000

IFU / REF MGES 1 / REV A / 2025-07-22

Note: Significant changes are indicated by dotted lines in the margin.
A change history can be found at the end of the instructions for use.

Milenia Biotec GmbH

Versailler Str. 1

35394 Gießen

Germany

Tel.: +49 641 94 8883-0

Fax: +49 641 94 8883-80

E-mail: info@milenia-biotec.de

Web: www.milenia-biotec.com

Table of Contents

English	8
Explanation of Symbols	10
Warnings and Precautions.....	10
Materials Supplied, Storage and Stability.....	11
Materials Required but not Supplied	11
Intended Use	11
Method	12
Test Procedure.....	12
Additional Products Available.....	13
Contact.....	13
Change History	14

Explanation of Symbols


Symbol	Explanation	Symbol	Explanation
	<i>Manufacturer</i>		Contains sufficient for <n> tests
	<i>Catalogue number</i>		Consult <i>instructions for use</i> or consult electronic <i>instructions for use</i>
	<i>Batch code</i>		Unique device identifier

Warnings and Precautions

- The product is intended for single use.
- Protect the materials from possible contamination. The primary packaging must be closed when not in use.
- Wear gloves during the work steps to reduce the risk of contamination.
- Waste disposal must be carried out in accordance with local regulations.
- For trained personnel.

Materials Supplied, Storage and Stability

Materials Supplied

REF	Content	Description	No. of Tests
MGES 1	2 x 500 Reaction Tube 1,5 mL 1 x 1000 Inner Reaction Tube	<u>Material</u> : 100% PP, transparent <u>Material</u> : 100% PP, transparent	1000

Storage and Stability

Components	Preparation	Storage	Shelf Life
Reaction Tube 1,5 mL	-	room temperature	-
Inner Reaction Tube	-	room temperature	-

Materials Required but not Supplied

In addition to the reaction tubes supplied with the kit, the following reagents and items are required:

- Extraction buffer MGSDDB
- Centrifuge (> 10.000g)
- Pipettes (100 µl)
- Pipette tips (containing protective filters)


Intended Use

The Milenia GenLine Extraction System is a universal kit designed for the extraction of analytes from swabs.

Method

The Milenia GenLine Extraction System consists of an inner and an outer tube. Swabs with sample material are placed in the inner vessel. The sample material is separated from the swab by adding a buffer. The solution with sample material is separated from the swab by centrifugation and collected in the outer tube.



Test Procedure

Caution: When taking samples, make sure that the swab is not contaminated with too much biomass (biofilm/mucus) or oily soiling. This can lead to inhibition of the PCR!

Sample processing :

1. Place the reaction vessel (A) in the reaction vessel (B).
2. Immediately after taking the sample, transfer the swab to the reaction vessel (A) and rotate it vigorously along the inner wall of the vessel at least 10 times.
3. Cut off the shaft and add 100 μ l swab detection buffer (MGSDB) to the swab. Allow the buffer to run slowly over the remaining shaft onto the swab. Close the reaction vessel (A) and incubate for at least 2 min at RT.
4. Centrifuge the samples for 1 min ($> 10,000g$). Position the lid against the direction of rotation. Discard the reaction tube (A) and close the reaction tube (B).

5. Resuspend the eluate for approx. 1 min using a vortexer. Finally, remove 2 μ l and use for the respective PCR detection!

Note: The eluate in the reaction vessel (B) can be stored at -20°C and analyzed at a later time.

Additional Products Available

REF	Produktname	Testanzahl	Beschreibung
MGUP 1	PCR Universal Module	50 tests	Lateral flow dipsticks for the detection of Genomic Amplificates (beer germs).
MGScLP 1	Lactobacillus/ Pediococcus-Screen	48 tests	Reagents for the detection of beer spoiling <i>Lactobacillus sp.</i> and <i>Pediococcus sp.</i> .
MGScMP 1	Megasphaera/ Pectinatus-Screen	48 tests	Reagents for the detection of beer spoiling <i>Megasphaera sp.</i> and <i>Pectinatus sp.</i> .
MGScHOR 1	Hop Resistance Gene-Screen	48 tests	Reagents for the detection of hop resistance genes.
MGSDb	Extraction Buffer	100 reactions	Buffer for purification of swab samples with MGES 1.

Contact

For additional information and support, please contact us:

Milenia Biotec GmbH

Versailler Str. 1

35394 Gießen

Germany

Tel.: +49 641 94 8883-0

Fax: +49 641 94 8883-80

E-mail: info@milenia-biotec.de

Web: www.milenia-biotec.com

Change History

Date	Revision	Cause of Revision
2025-07-22	IFU/REF MGES 1/ REV A /2025-07-22	Creation of the IFU