# Milenia GenLine

## ***Megasphaera/Pectinatus-Screen***

---

Deutsch

---

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DE: Seite 1-13

Reagents for the detection of beer spoiling  
*Megasphaera sp.* and *Pectinatus sp.*

EN: Page 14-26

---

**REF** MGScMP 1\*  48

---

IFU / REF MGScMP 1 / REV D / 2025-07-10

**\* Attention:**

This kit contains only reagents used for amplification.

For the **detection** of amplicates the following additional kit is required:

**Milenia GenLine PCR-Universal module REF MGUP 1 (50 tests)**

Note: Significant changes are indicated by dotted lines in the margin.  
A change history can be found at the end of the manual.

**Milenia Biotec GmbH**

Versailler Str. 1

35394 Gießen

Germany

Tel.: +49 641 94 8883-0

Fax: +49 641 94 8883-80

E-mail: [info@milenia-biotec.de](mailto:info@milenia-biotec.de)

Web: [www.milenia-biotec.com](http://www.milenia-biotec.com)

## Table of Contents


|                                                                 |    |
|-----------------------------------------------------------------|----|
| English .....                                                   | 14 |
| Explanation of Symbols .....                                    | 16 |
| Warnings and Precautions.....                                   | 16 |
| Materials Supplied, Storage and Stability.....                  | 17 |
| Safety Data Sheets (SDS) .....                                  | 18 |
| Materials Required but not Supplied .....                       | 18 |
| Intended Use .....                                              | 19 |
| Method and Test Principle.....                                  | 19 |
| Procedure of the Test – A short Overview.....                   | 19 |
| Preparation of Samples and PCR Reagents.....                    | 20 |
| a) Preparation of PCR Reaction Mix.....                         | 20 |
| b) Sample preparation from liquid culture media .....           | 21 |
| c) Sample preparation of single colonies from solid media ..... | 21 |
| d) Sample preparation of swab samples .....                     | 21 |
| Test Procedure.....                                             | 22 |
| a) PCR Reaction.....                                            | 22 |
| b) Detection of PCR Products .....                              | 22 |
| Interpretation of Results.....                                  | 23 |
| a) Positive Result .....                                        | 23 |
| b) Negative Result .....                                        | 23 |
| c) No Analysis Possible .....                                   | 24 |
| Troubleshooting .....                                           | 24 |
| References .....                                                | 24 |
| Additional Products Available.....                              | 25 |
| Contact.....                                                    | 25 |
| Change History .....                                            | 26 |

## Explanation of Symbols

| Symbol                                                                            | Explanation             | Symbol                                                                            | Explanation                                                                           |
|-----------------------------------------------------------------------------------|-------------------------|-----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  | <i>Manufacturer</i>     |  | Contains sufficient for <n> tests                                                     |
|  | Use-by date             |  | Consult <i>instructions for use</i> or consult electronic <i>instructions for use</i> |
|  | <i>Batch code</i>       |  | Unique device identifier                                                              |
|  | <i>Catalogue number</i> |  | Temperature limit                                                                     |

## Warnings and Precautions

- Read this instruction for use carefully prior running the test. Follow the sequence of the work steps.
- Pay attention to the expiration dates of the individual components of the test.
- Complete thawing of all components prior use is required.
- Do not exchange components with different Lot numbers.
- Store all reagents at **-15 °C to -25 °C** in the original packaging.
- **Immediately after use**, store all components at **-15 °C to -25 °C**.
- The disposal of waste materials must be carried out according to current local regulations.
- Wear laboratory coat and protective gloves.
- In order to avoid contaminations, run the different steps of the test at different locations.
- The place where the reaction mixes are prepared must be separated from the place where cultures are handled and where the evaluation of PCR is done.
- Decontamination of the work area and the equipment used must be carried out regularly!

**Figure 1:** Separation of the testing procedure in different working areas

## Materials Supplied, Storage and Stability

### Materials Supplied

| REF      | Content                       | Part.-No. | No. of Tests |
|----------|-------------------------------|-----------|--------------|
| MGScMP 1 | 2 x PCR-Polymerase Master Mix | MGPOL     | 48           |
|          | 2 x PCR-Primer Mix            | MGPMMP    |              |
|          | 1 x Positive Kit-Control      | MGMPPC    |              |
|          | 1 x PCR-Water *               | MGH2O     |              |

\* The PCR-Water (MGH2O) is used as diluent and negative control.

### Storage and Stability

| Kit Component                     | Content | Preparation  | Storage           | Shelf Life        |
|-----------------------------------|---------|--------------|-------------------|-------------------|
| PCR-Polymerase Master Mix (MGPOL) | 300 µL  | ready-to-use | -15 °C bis -25 °C | until expiry date |
| PCR-Primer Mix (MGPMMP)           | 280 µL  | ready-to-use | -15 °C bis -25 °C | until expiry date |
| Positive Kit Control (MGMPPC)     | 250 µL  | ready-to-use | -15 °C bis -25 °C | until expiry date |
| PCR-Water (MGH2O)                 | 1000 µL | ready-to-use | -15 °C bis -25 °C | until expiry date |

## Safety Data Sheets (SDS)

Safety Data Sheets (SDS) are available on request and alternatively under the respective product on our website ([www.milenia-biotec.com](http://www.milenia-biotec.com)).

[SDS REF MGScMP 1 - \*Megasphaera/Pectinatus\*-Screen](#)

## Materials Required but not Supplied

In addition to the components supplied with the kit, the following reagents and items are required:

- **Milenia GenLine PCR-Universal module (MGUP 1)**. Contains test strips and chase buffer.
- **Thermocycler Labcycler 48s, (MLCY)**, available at Milenia Biotec GmbH
- Deep Freezer: -20 °C
- Pipets 2- 20 µL, 20- 200 µL, 100- 1,000 µL
- Sterile (DNA-free) pipet tips (with contamination protection for PCR)
- Sterile (DNA-free) reaction tubes 1.5 mL
- Sterile (DNA-free) PCR-tubes 0.2 mL; recommended: Reaction vessels from Sorensen Bioscience Inc. due to thin wall thickness for optimum temperature transfer
- 96-well microtiter plate
- (Desktop) centrifuge for reaction- and PCR-tubes
- Vortex mixer
- Tube holder for reaction- and PCR-tubes, e.g.: isofreeze PCR tube rack

The following reagents and aids are also required **for hygiene monitoring with swabs**:

- **Swab Detection Buffer (REF MGSDB)**
- **Milenia GenLine Extraction System (REF MGES 1)**
- Centrifuge (> 10,000g)
- Swab

## Intended Use

The Milenia GenLine *Megasphaera/Pectinatus*-Screen is designed for the qualitative detection of beer spoiling *Megasphaera sp.* and *Pectinatus sp.*.

## Method and Test Principle

Milenia GenLine *Megasphaera/Pectinatus*-Screen is a rapid detection system based on PCR (polymerase-chain-reaction). The screening test is combined with the evaluation kit MGUP 1 which allows the user a optical evaluation of the PCR reaction via lateral flow test strips.

In a first part of the test, labelled specific primers (short, single stranded DNAs) are used to amplify specific DNA fragments of *Megasphaera sp.* and *Pectinatus sp.*. Thus, the resulting PCR products carry the labels of the incorporated primers. In addition to the target gene, a control gene, which is also present in the PCR mixes, is amplified to make sure that the PCR process works properly.

In a second part of the test, the created PCR products are detected by a lateral flow test strip. The PCR products bind via their labels to antibodies attached to the solid phase of the strip and to colloidal gold particles which are coated with an antibody directed against the other label. Thus a "molecular sandwich" is formed and becomes visible as a line on the test strip.

## Procedure of the Test – A short Overview

The entire procedure of the test can be divided into 3 different steps.

1. Preparation of PCR reagents and sample
2. PCR reaction
3. Evaluation by a simple lateral flow test strip

**Figure 2:** Presentation of how the procedure is carried out

## Preparation of Samples and PCR Reagents

The preparation of samples and reagents includes the preparation of the PCR reaction mix and the samples to be analysed.

### a) Preparation of PCR Reaction Mix

The precise and contamination free preparation of the mix is the prerequisite of a proper function of the test. It is recommended to label all reaction tubes prior pipetting the samples. In case a positive (MGMPPC) and a negative (MGH2O) control is included in the run, add 2 additional tubes to the number of samples introduced in the test run.

Procedure:

1. According to the number of samples (plus positive and negative control) pipette **PCR Polymerase Master Mix** into a sterile 1.5 mL reaction tube (see Table 1).
2. Add an appropriate volume of **PCR Primer Mix** (see Table 1).
3. **Mix** the two solutions (reaction mix) **well** (pipet, Vortex).
4. Prepare **sterile 0.2 mL PCR reaction tubes** according to the number of planned reactions.
5. Pipet **23 µL** of the prepared reaction mix into each PCR reaction tube.
6. **Immediately close** the PCR-reaction tubes.  
The sample is added at a separate workstation.

**Table 1:** Overview over the volumes to be pipetted depending on the number of PCR reaction tubes.

| Number of<br>PCR Reaction Tubes | Pipette volume<br>PCR Polymerase Master Mix | Pipette volume<br>PCR Primer Mix |
|---------------------------------|---------------------------------------------|----------------------------------|
| 1                               | 13 µL                                       | 11 µL                            |
| 5                               | 65 µL                                       | 55 µL                            |
| 10                              | 130 µL                                      | 110 µL                           |
| 20                              | 260 µL                                      | 220 µL                           |

## b) Sample preparation from liquid culture media

No sample preparation is needed if the analysis is done from liquid culture media.

The following media based on MRS can be used for the test:

|                |                     |                           |
|----------------|---------------------|---------------------------|
| <u>Medium:</u> | NBB-B (Bouillon)    | – Döhler GmbH (Darmstadt) |
|                | NBB-B-AM (Bouillon) | – Döhler GmbH (Darmstadt) |
|                | NBB-C (Concentrate) | – Döhler GmbH (Darmstadt) |

2 µL from the enriched culture are transferred to the prepared PCR reaction mix.

### **Note:**

Very turbid enriched cultures contain large amounts of biological material. This can have a negative impact on the test or even may inhibit it completely. Therefore, such samples must be pre-diluted 1:10 or 1:100 prior testing.

## c) Sample preparation of single colonies from solid media

### Procedure:

1. Touch a colony on the agar plate with a pipet tip. Do not wipe across the colony.
2. Resuspend the cells in 100 µL sterile water.
3. Add 2 µL of the cell suspension to the PCR reaction mix.

### **Note:**

If too many cells are added to the mix, this may have a negative impact on the test. Therefore, a careful light touch in step 1 is recommended.

## d) Sample preparation of swab samples

### Procedure:

1. Sampling with a swab at process-critical points.
2. Detach the sample material from the swab (e.g. using the Milenia GenLine extraction system).
3. 2 µL of the detached sample material is added to the prepared PCR reaction mix.

We recommend the Milenia GenLine Extraction System (REF: MGES 1) and the Swab Detection Buffer (REF: MGSDDB) for purification of the swab samples. Further information can be found on the [Milenia GenLine Extraction System](#) product page.

## Test Procedure

The test procedure describes the start and termination of the PCR reaction as well as the subsequent lateral flow based detection of PCR products by a universal Test Strip (Milenia GenLine PCR Universal Module MGUP 1).

### a) PCR Reaction

1. After addition of the sample to the PCR tube containing the reaction mix, close the tube and mix well.
2. Place all PCR tubes into a minicentrifuge and centrifuge the samples briefly.
3. All liquid should be free of air bubbles and located on the bottom of the PCR tubes.
4. Place all PCR tubes into a thermocycler.
5. Close the lid of the thermocycler.
6. Start the species-specific program "MP-Screen" within the folder „Bier“.
7. After completion of the program, stop the program and remove the PCR tubes from the thermocycler.

When the program is finished, the thermocycler holds a temperature of 8 °C. PCR products must be stored in a refrigerator (app. 4 °C) and should be analyzed subsequently (max. 24 h later).

### b) Detection of PCR Products

For the detection of PCR products, the **Milenia GenLine PCR-Universal Module REF MGUP 1** is required.

It is strongly recommended to use pipet tips with filter and separate pipets for the addition of the sample to be analyzed and the evaluation of PCR products to minimize the risk of DNA contamination. Regular decontamination of the pipettes and the separation of the different working steps at different locations are recommended.

#### Procedure:

1. According to the number of tests, pipet 80 µL chase buffer (MGCBA) each into clean wells of a 96-well microtiter plate.
2. Label each test strip (MGDS2A) needed for the analysis and take it aside.
3. Pipet 2 µL of each PCR product onto the sample application pad of a test strip.


4. Place the test strip with the sample application pad down into the chase buffer of the wells of the microtiter plate and incubate for 5 min at room temperature. Read the results visually after exactly 5 min.
5. A documentation by taking an image of the strips with a camera is recommended.

**Note:**

Results must be read exactly after 5 minutes! After prolonged time unspecific, weak signals may occur on the test strip which are not associated with positive results.

## Interpretation of Results

Tests are interpreted visually. Potential test results are shown in the figure below.

**Figure 3:** Potential test results for visual interpretation

The control line [C] confirms the proper function of the test strip. The test line 1 [T1] is associated with the presence or absence of the target gene. The test line 2 [T2] is linked with the proper function of the PCR.

### a) Positive Result

The prerequisite for a valid positive result is the appearance of a distinct test line [T1] (test strip **1 and 2**). If high levels of the target gene are introduced into the test, the amplification control can vanish. Nevertheless, there will be a strong signal on test line T1. Therefore, no line may be visible at the test line [T2]. Such results should be read positive (test strip **2**).

### b) Negative Result

A valid negative result requires a distinct visible test line [T2]. This confirms that the PCR has worked properly. The test line [T1] is not visible (test strip **3**).

**Note:**

If the line of the PCR control [T2] is weakly visible, an inhibition of the PCR is very likely (test strip **3\***). In this case it is recommended to rerun the samples after predilution.

### c) No Analysis Possible

If neither the test line [T1] nor the PCR control [T2] is visible, no result can be reported (test strip **4**). In case only the line at the location of the control line [C] is visible, either the PCR product has not been transferred successfully to the strip or the amplification control has not performed properly.

If **no line** is visible on the test strip (test strip **5**), please contact the customer service at Milenia Biotec GmbH.

## Troubleshooting

### The negative control gives a positive result.

**Contaminations are a high risk when working with nucleic acids. It is very difficult to find out the reasons for contamination of materials of cellular origin or of purified DNA. Therefore, the following preventive actions are strongly recommended:**

- Check that the chosen PCR protocol of the thermocycler was correct.
- Decontaminate the workspace where the preparation of the reaction mix is done.
- Decontaminate the pipets.
- Use unopened, sterile plastic materials (pipet tips/ 1.5 mL reaction vials/PCR vials).

After decontamination, a test run should be performed in order to check the effectiveness of the preventive actions. It is recommended to run a positive and a negative control. The negative control **must** show a negative result.

If the negative control shows a positive result, this may be due to contamination of the reagents. In this case new reagents (PCR-Water, PCR-Polymerase Master Mix (2x) and PCR-Primer Mix (2x)) must be used. In the subsequent control run the negative control (PCR-Water, MGH2O) must be negative.

If the negative result shows a positive result again, please contact customer service at Milenia Biotec GmbH.

## References

1. Koob, J et al., Molekularbiologischer Nachweis bierschädigender Organismen, Der Weihenstephaner, October 2015
2. Koob, J: „A New Approach to Bacterial Identification in a Point-of-Care Format“; European Brewery Convention. Mundiconvenius, Porto, Portugal (2015)

## Additional Products Available

| REF       | Produktname                                      | Testanzahl     | Beschreibung                                                                                    |
|-----------|--------------------------------------------------|----------------|-------------------------------------------------------------------------------------------------|
| MGUP 1    | <a href="#">PCR Universal Module</a>             | 50 tests       | Lateral flow dipsticks for the detection of Genomic Amplificates (beer germs).                  |
| MGScHOR 1 | <a href="#">Hop Resistance Gene-Screen</a>       | 48 tests       | Reagents for the detection of hop resistance genes.                                             |
| MGScLP 1  | <a href="#">Lactobacillus/Pediococcus-Screen</a> | 48 tests       | Reagents for the detection of beer spoiling <i>Lactobacillus sp.</i> and <i>Pediococcus sp.</i> |
| MGES 1    | <a href="#">Extraction System</a>                | 1000 reactions | Universal kit for the isolation of analytes from swabs for further processing.                  |
| MGSDb     | <a href="#">Extraction Buffer</a>                | 100 reactions  | Buffer for purification of swab samples with MGES 1.                                            |

## Contact

For additional information and support, please contact us:

**Milenia Biotec GmbH**

Versailler Str. 1

35394 Gießen

Germany

Tel.: +49 641 94 8883-0

Fax: +49 641 94 8883-80

E-mail: [info@milenia-biotec.de](mailto:info@milenia-biotec.de)

Web: [www.milenia-biotec.com](http://www.milenia-biotec.com)

## Change History

| Date       | Revision                                      | Cause of Revision                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|------------|-----------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 2017-03-01 | MGScHOR 1/ <b>B</b> /2017-03-01               | Insertion of english description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| 2019-06-21 | MGScLP 1/ <b>C</b> /2019-06-21                | Adjustment of the MGUP 1 quantity                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 2025-07-10 | IFU/REF MGScMP 1/REV <b>D</b> /<br>2025-07-10 | IFU created in a new design and revised in the following chapters: <ul style="list-style-type: none"><li>- Explanation of Symbols (<b>UDI</b> symbol and symbol <b>Contains sufficient for &lt;n&gt; tests</b>)</li><li>- Materials Supplied, Storage and Stability</li><li>- Safety Data Sheets (new chapter created)</li><li>- Method and Test Principle</li><li>- Preparation of Samples and PCR Reagents (point <b>d</b>) <b>Sample preparation of swab samples</b> added)</li><li>- Test Procedure (note added)</li><li>- Additional Products available (new chapter created)</li></ul> |