Milenia QuickLine HIT Control

English

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Milenia QuickLine HIT Control **REF** MQCOHIT 1  20

IFU / REF MQCOHIT 1 / REV E / 2025-06-04

Note: Significant changes are indicated by dotted lines in the margin.
A change history can be found at the end of the manual.

Milenia Biotec GmbH
Versailler Str. 1
35394 Gießen
Germany
Tel.: +49 641 94 8883-0
Fax: +49 641 94 8883-80
E-mail: info@milenia-biotec.de
Web: www.milenia-biotec.com

Table of Contents

English

Explanation of Symbols	3
Warnings and Precautions.....	3
Intended Use	4
Materials Supplied, Storage and Stability.....	4
Additional Materials Required (not included).....	4
Volume Required.....	4
Preparation and Procedure.....	5
Interpretation of Results.....	5
Internal Quality Control Process	6
Additional Products available	6
Contact.....	6
Change History	7

Explanation of Symbols

Symbol	Explanation	Symbol	Explanation
	Manufacturer		<i>In vitro</i> diagnostic medical device
	Use-by date		Consult instructions for use or consult electronic instructions for use
	Batch Code		Unique Device Identification
	Catalogue number		Control
	Contains sufficient for <n> tests		Positive Control
	Temperature limit		Negative Control
	CE Marking		Contains biological material of animal origin

Warnings and Precautions

- Read this instruction manual and the instruction manual of Milenia QuickLine HIT-Test (MQHIT 1) carefully before performing the assay.
- Follow the sequence of the work steps.
- Aliquoted control should be thawed completely before usage.
- Allow all reagents reach room temperature before use.
- Handle Milenia QuickLine HIT Control reagents carefully.
- Wear lab coat and protective gloves.
- The disposal of waste must be carried out according to local regulations.
- Animal materials which are used in the reagent are free from infection pathogens and are derived from animals kept in Germany.
- Human materials which are used in the reagent were tested and found negative for: HIV 1/2, HbsAg, HCV, HIV-1 (NAT) and RPR by FDA approved methods.
- For professional users.

Safety data sheets are available on request.

Intended Use

The Milenia QuickLine HIT Control is used for internal quality control of the corresponding Milenia QuickLine HIT-Test. The Milenia QuickLine HIT Control kit contains a negative and a positive sample. It is used to determine test validity. The product is for in vitro diagnostic use.

Materials Supplied, Storage and Stability

Component	Cat.-No.	Content	Preparation	Storage (unopened)	Shelf Life
QuickLine HIT Positive Control	MQHITC1	1 vial with lyophilisate (50 tests)	lyophilisate; reconstitute with 0,25 mL distilled water not less than 30 min before usage	2-8°C	lyophilised at 2°C to 8°C stable until the imprinted expiry date
QuickLine HIT Negative Control	MQHITC2	1 vial with lyophilisate (50 tests)	lyophilisate; reconstitute with 0,25 mL distilled water not less than 30 min before usage	2-8°C	lyophilised at 2°C to 8°C stable until the imprinted expiry date

Storage and component stability after reconstitution	Storage	Shelf Life
Control reconstituted in aqua dest.	2-8°C	1 month
Control is aliquoted immediately after reconstitution	-20°C or below	3 months

Additional Materials Required (not included)

- Pipette and tips: 250 µL
- Distilled water
- Milenia QuickLine HIT-Test (MQHIT 1)
- Multi-axe rotating mixer

Volume Required

5 µL Milenia QuickLine HIT Positive/Negative Control is required per test unit (MQHIT 1).

Preparation and Procedure

a) Please follow these instructions to prepare the control material for first-time use

- Bring lyophilised Milenia QuickLine HIT Positive and Negative Control to room temperature (18°C to 25°C) before dissolving
- **Add 250 µL distilled water** to each vial
- Close control vials
- Mix the Milenia QuickLine HIT Positive and Negative Control thoroughly with e.g. a Vortex Mixer. Afterwards put the control vials for at least 30 minutes on a multi-axe rotating mixer (if possible).
- The Milenia QuickLine HIT Control is now ready-to-use.
- For long-term storage (max. 3 month), it is recommended to divide, and aliquot Milenia QuickLine HIT Control material and freeze at - 20 °C or below until it is needed. Do not refreeze.

b) Use of Milenia QuickLine HIT Control:

- **Add 5 µL** Milenia QuickLine HIT Positive or Negative Control in the MQHIT testunit. Before use, please consult instructions for use of Milenia QuickLine HIT (Ref MQHIT 1).

Please note: The control material should be brought to room temperature, if it was stored frozen or cooled (at least 30 min on a multi-axe rotating mixer) before using it in the MQHIT Test.

Interpretation of Results

Fig. 1 Positive Control

- Milenia QuickLine HIT **Positive** Control is considered positive, if the test line (T) shows the **same** or a **darker** shade of color than the given line on the Evaluation Card.


Fig. 2 Negative Control

- Milenia QuickLine HIT **Negative** Control is evaluated negative if the test line (T) does not appear or the line shows a **weaker shade of color** than the given line on the Evaluation Card.

Internal Quality Control Process

Each laboratory should establish quality control procedures that are tailored to their individual requirements and are in compliant with relevant regulations or the terms of approval, including corrective actions which must be taken if the results are outside of target ranges.

We recommend running a positive and a negative control when starting to use a new Lot.

In addition, the specifications written in the instructions of Milenia QuickLine HIT test (REF MQHIT 1) is subject to the evaluation of Milenia Quickline HIT Control.

Additional Products available

Product Name	Order No.	Content	Description
Milenia QuickLine HIT	MQHIT 1	20 tests	Qualitative detection of IgG antibodies against PF4/polyanion-complexes in human citrated plasma or serum.

Contact

For additional information and support, please contact us:

**Milenia Biotec GmbH**

Versailler Str. 1

35394 Gießen

Germany

Tel.: +49 641 94 8883-0

Fax: +49 641 94 8883-80

E-mail: info@milenia-biotec.de

Web: www.milenia-biotec.com

Change History

Date	Revision	Cause of Revision
2017-02-28	MQCOHIT 1 Rev B 2017-03-01	Registered Tradename changed from „Milenia Quickline HIT Positive Kontrolle“ to „Milenia Quickline HIT Positive Control“
2019-06-04	MQHIT 1_Z Rev C 2019-06-04	Addition of the Negative Control
2022-04-12	MQHIT 1_Z Rev D 2022-04-12	Adaption of symbols (cover and explanation)
2025-06-04	IFU / REF MQCOHIT 1 / REV E / 2025-06-04	• Symbols adapted in chapter <i>Explanation of Symbols</i>. • Design updated

Milenia QuickLine HIT Control

Deutsch

Qualitative HIT Control, exclusively designed for the
application in Milenia QuickLine HIT

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