Milenia QuickLine HIT

English

**Lateral Flow Immunoassay to detect IgG antibodies
against PF4/polyanion-complexes**

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 20

IFU / REF MQHIT 1 / REV H / 2025-06-04

Note: Significant changes are indicated by dotted lines in the margin.
A change history can be found at the end of the manual.

Milenia Biotec GmbH

Versailler Str. 1

35394 Gießen

Germany

Tel.: +49 641 94 8883-0

Fax: +49 641 94 8883-80

E-mail: info@milenia-biotec.de

Web: www.milenia-biotec.com


Table of Contents

English

Explanation of Symbols	3
Warnings and Precautions.....	3
Materials Supplied, Storage and Stability.....	4
Intended Use	4
Materials Required but not Supplied	4
Specimen Collection and Preparation	4
Volume Required.....	5
Summary and Explanation.....	5
Method and Test Principle.....	5
References	6
Test Performance.....	7
Interpretation of Results.....	8
Function Control	9
Assay Characteristics	9
Additional Products available	9
Contact.....	10
Change History	10

Explanation of Symbols

Symbol	Explanation	Symbol	Explanation
	Manufacturer		Do not re-use
	Use-by date		Consult instructions for use or consult electronic instructions for use
	Batch Code		Unique Device Identification
	Catalogue number		<i>In vitro</i> diagnostic medical device
	Do not use if package is damaged and consult instructions for use		Contains sufficient for <n> tests
	Temperature limit		CE Marking
	Contains biological material of animal origin		Contains biological material of human origin

Warnings and Precautions

- Read this instructions for use carefully before performing the test. Follow them step by step.
- Store reagents at 2 - 8° C.
- Do not use Test Unit if foil pouch is damaged.
- Protect Test Units from humidity.
- Do not interchange components of packages of different batch codes.
- All samples, used Test Units, pipettes, tips, and other contaminated material should be handled carefully and as potentially infectious.
- The disposal of waste materials must be carried out according to current local regulations.
- For *in vitro* diagnostic use only. These reagents are to be used only by certified medical laboratory personnel authorized by the laboratory.
- The reagents provided in these kits contain materials of human and/or animal origin. Whenever human blood product is required for the preparation of these reagents, approved methods are used to test the blood for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method

can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with safety precautions in the manipulation of these biological materials as if they were infectious.

- For professional users.

Materials Supplied, Storage and Stability

Components	Cat.-No.	Content	Preparation	Storage at	Shelf Life
HIT Test Units Membrane coated with polyclonal (goat) anti-human IgG antibodies and anti-ligand antibodies conjugated to gold particles	MQSHIT	20	ready to use	2-8°C Use immediately after opening the foil pouch!	in the closed foil pouch until expiry date
Sample Buffer Buffer, pH 7.4	MQBHIT	1 vial à 2 mL (white)	ready to use	2-8°C	until expiry date
Evaluation Cards	MQAHIT	20			

Material Safety Data Sheets are available on request (see also www.milenia-biotec.de).

Intended Use

- For *in vitro* diagnostic use.
- Milenia QuickLine HIT is designed for qualitative detection of IgG antibodies against PF4/polyanion-complexes in human citrated plasma or serum.

Materials Required but not Supplied

- Pipet and tips: 5 µL

Specimen Collection and Preparation

Serum or citrated plasma

Sample collection must be in conformity with the recommendations for haemostasis tests.

Blood (9 vol.) is collected in 0.109 M (i.e., 3.2 %) trisodium citrate anticoagulant (1 vol.) or in a tube without anticoagulant.

Centrifugation: Within two hours after the sample collection

Do not use samples with visible fibrin strands. Such samples need to be centrifuged for a second time until no fibrin strands are visible.

Volume Required

- Citrated plasma or serum : 5 µl

Summary and Explanation

Heparin induced thrombocytopenia (HIT) is a life-threatening disease associated with exposure to unfractionated or (less commonly) low-molecular-weight heparin, drugs which are used very frequently for various indications. HIT II occurs in up to 5 % of patients on heparin. This prevalence varies with patient population, the nature of the disease and the type of heparin, thus the anticoagulant. Cardiac and major orthopedic surgeries expose a higher risk than other clinical situations when heparin is used.

HIT is caused by antibodies directed against a complex of platelet factor 4 (PF4) and polyanions such as heparin. The pathogenic and thrombogenic species is immunoglobulin IgG. There is evidence for a correlation between antibody concentration and clinical severity. Platelet activation and massively enhanced thrombin generation promote venous and/ or arterial thromboembolism. The rate of severe patient injury (e.g. amputations) and mortality is very high.

Patients with HIT II require alternative anticoagulation, but the drugs used in this setting have a high rate of hemorrhagic complications. There is no antidote available. Monitoring with standard methods is either not reliable or impossible. On top of this the drug costs are significantly higher than for heparins. Alternative anticoagulants should be reserved for patients with high clinical probability for HIT and with confirmed antibodies. Excluding HIT-antibodies is an important step towards selecting the most appropriate anticoagulation strategy.

Method and Test Principle

The Milenia QuickLine HIT is a lateral flow immunoassay designed for qualitative determination of antibodies against PF4/polyanion-complexes in human citrated plasma or serum.

The sample is pipetted in the sample application port of the test unit and the sample buffer is added immediately. The buffer forces the sample to migrate through the membrane of the test unit.

Patient's IgG-antibodies are captured by the PF4/polyanion-complex and bound to immobilized anti-human IgG on the nitrocellulose membrane. The complex is detected by intensely colored gold nanoparticles.

The presence of PF4/polyanion-complexes becomes visible at a colored test line (T). This test line - if it is positive - has to be the same or a darker shade of color than the given one on the provided evaluation card.

The surplus of gold particles continues to migrate through the membrane and is captured at the control line (C) by specific antibodies. A well visible line appears after incubation time.


References

1. Greinacher A. Heparin-induced thrombocytopenia. J Thromb Haemost. 2009; 7 Suppl 1:9-12. Review.
2. Wilke T, Tesch S, Scholz A, Kohlmann T, Greinacher A. The costs of heparin-induced thrombocytopenia: a patient-based cost of illness analysis. J Thromb Haemost. 2009;7(5):766-73.
3. Bakchoul T, Giptner A, Najaoui A, Bein G, Santoso S, Sachs UJ. Prospective evaluation of PF4/heparin immunoassays for the diagnosis of heparin-induced thrombocytopenia. J Thromb Haemost. 2009;7 (8):1260-5.
4. Sachs UJ, von Hesberg J, Santoso S, Bein G, Bakchoul T. Evaluation of a new nanoparticle-based lateral-flow immunoassay for the exclusion of heparin-induced thrombocytopenia (HIT). J Thromb Haemost. 2011; 106(6).
5. Warkentin TE, Greinacher A, Koster A, Lincoff AM. American College of Chest Physicians. Treatment and prevention of heparin-induced thrombocytopenia: American College of Chest Physicians Evidence-Based Clinical Practice Guidelines (8th Edition). Chest. 2008 un;133(6 Suppl):340S-380S.

Test Performance

Important Notes:

- Before opening bring the test unit and sample buffer to room temperature (18 – 28 °C)!
- Use the test unit immediately after opening. Do not store at room temperature!
- Do not use the test unit without patient sample, this can lead to crossreactions with the buffer.

Step		
1. Sample	Add 5 µL serum or citrated plasma to the round sample port of the test unit.	
2. Buffer	Add 2 drops Sample Buffer Note: Do not contaminate the vial with sample material! Avoid air bubbles!	
3. Incubation	10 minutes at room temperature (18 - 28°C)	
4. Evaluation	The result can only be interpreted visually by the provided Evaluation Card: The test is interpreted positive if the test line shows the same or a darker shade of color after incubation time than the given line in the evaluation card (s. Results Interpretation and Functional Control) The result must be interpreted immediately after incubation time.	

Interpretation of Results

Hold the colored line of the Evaluation Card (as shown below) **close beside** the test line of your test unit. (See figures 1, 2)



Figure 1

Figure 2

The **result** will be evaluated **negative**, if the **Test Line-T** shows a **lighter shade of colour** than the given line on the Evaluation Card. If no **T-Line** is visible the result is also **negative**.

The **result** will be evaluated **positive**, if the **Test Line-T** shows the **same (Fig. 1) or a darker (Fig.2) shade of colour** than the given line on the Evaluation Card.

The **Control Line-C** must be visible. (see **Function Control** of the test unit) for a valid result.


Note: Results should be verified concerning the entire clinical status of the patient. Additionally every decision for therapy should be taken individually.

Function Control

In any case, the Control line (C) has to appear!

This line is only used as functional control and is no result.

If the control line is not detectable after 10 minutes incubation time, the result is invalid and the analysis has to be repeated with a new test unit!

For the quality assurance, it is recommended to run a control (MQCOHIT 1) every time a new Lot is used. Additional controls should be done regular according to the medical necessity and the frequency of examination of patient samples.

Assay Characteristics

Sample Material:	5 µL serum or citrated plasma
Incubation Time:	10 minutes at room temperature (18 – 28 °C)
Normal Range	negative
Diagnostic Sensitivity:	100%
Diagnostic Specificity:	93%
Interference Factors:	Lipemia, bilirubin and haemolysis do not show any interference influence on the results. Biotin shows no interference up to the maximum concentration tested (1200 ng/mL).

Additional Products available

Product Name	Order No.	Content	Description
Milenia QuickLine HIT Control	MQCOHIT 1	2 x 50 tests	The positive and negative controls are designed for internal quality control of Milenia QuickLine HIT tests.

Contact

For additional information and support, please contact us:

**Milenia Biotec GmbH**

Versailler Str. 1

35394 Gießen

Germany

Tel.: +49 641 94 8883-0

Fax: +49 641 94 8883-80

E-mail: info@milenia-biotec.de

Web: www.milenia-biotec.com

Change History

Date	Revision	Cause of Revision
2013-07-15	MQHIT 1_Z Rev E 2013-07-15	Plasma can be used as well as sample material.
2014-02-10	MQHIT 1_Z Rev F 2014-02-10	Plasma specified as citrated plasma
2020-07-15	MQHIT 1_Z Rev G 2020-07-15	Addition of biotin to Interference Factors
2025-06-04	IFU / REF MQHIT 1 / REV H / 2025-06-04	• Symbols adapted in chapter <i>Explanation of Symbols</i>. • Design updated

Milenia QuickLine HIT

Deutsch

Lateral Flow Immunoassay to detect IgG antibodies
against PF4/polyanion-complexes

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Referenzen

1. Greinacher A. Heparin-induced thrombocytopenia. J Thromb Haemost. 2009; 7 Suppl 1:9-12. Review.
2. Wilke T, Tesch S, Scholz A, Kohlmann T, Greinacher A. The costs of heparin-induced thrombocytopenia: a patient-based cost of illness analysis. J Thromb Haemost. 2009; 7(5):766-73
3. Bakchoul T, Giptner A, Najaoui A, Bein G, Santoso S, Sachs UJ. Prospective evaluation of PF4/heparin immunoassays for the diagnosis of heparin-induced thrombocytopenia. J Thromb Haemost. 2009; 7 (8):1260-5.
4. Sachs UJ, von Hesberg J, Santoso S, Bein G, Bakchoul T. Evaluation of a new nanoparticle-based lateral-flow immunoassay for the exclusion of heparin-induced thrombocytopenia (HIT). J Thromb
5. Warkentin TE, Greinacher A, Koster A, Lincoff AM; American College of Chest Physicians. Treatment and prevention of heparin-induced thrombocytopenia: American College of Chest Physicians Evidence-Based Clinical Practice Guidelines (8th Edition). Chest. 2008 Jun;133(6 Suppl):340S-380S.

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