



Literature References citing the use of Milenia HybriDetect

Date: February 3rd, 2021

The Milenia HybriDetect is a universal lateral flow development platform that allows researchers all over the world to develop their own rapid test. Today the HybriDetect is frequently cited in more than 140 peer reviewed publications, including high impact journals such as Nature or Science.

Total number of publications citing the HybriDetect: 145

The HybriDetect test strips have developed into a valuable development tool, which is perfectly compatible with DNA amplification methods. The robust and easy-to-handle dipsticks are an interesting alternative for the detection of amplification products under low resource settings. Therefore, many researchers use the HybriDetect to underline the Point-of-Care compatibility of their unique assay. Nucleic acid amplification techniques can be multiplexed, highly specific and extremely sensitive. The combination of DNA amplification and Lateral Flow is a rapidly evolving field.

Over 85% of all papers citing HybriDetect use isothermal DNA amplification

So far the HybriDetect has successfully combined with the polymerase chain reaction (PCR), isothermal amplification methods like the loop mediated isothermal amplification (LAMP) or recombinase polymerase amplification (RPA). The 2020 pandemic, caused by SARS-CoV-2, brought molecular, Point-of-Care compatible diagnostic into public and scientific focus.

2020- the year of CRISPR/Cas-based detection strategies: 21 publications

Table of contents:

<i>Polymerase Chain Reaction (PCR)</i>	page 2
<i>Loop mediated isothermal Amplification (LAMP)</i>	page 3
<i>Recombinase Polymerase Chain Reaction (RPA)</i>	page 7
CRISPR/Cas-based Nuclease Assays.....	page 13
Special Applicationa (tHDA, CPA, RCA,).....	page 17



Polymerase Chain Reaction (PCR)

1. Loose, F. N., Breitbach, A., Bertalan, I., Ruster, D., Truyen, U., & Speck, S. (2020). Diagnostic validation of a rapid and field-applicable PCR-lateral flow test system for point-of-care detection of cyprinid herpesvirus 3 (CyHV-3). *PLoS ONE*, 15(10 October). <https://doi.org/10.1371/journal.pone.0241420> **multiplex**
2. Seo, H., Kil, E. J., Fadhila, C., Vo, T. T. B., Auh, C. K., Lee, T. K., & Lee, S. (2020). Rapid diagnosis of two marine viruses, red sea bream iridovirus and viral hemorrhagic septicemia virus by PCR combined with lateral flow assay. *VirusDisease*, 31(3), 251–256. <https://doi.org/10.1007/s13337-020-00577-z>
3. Wint, N. Y., Han, K. K., Yamprayooswat, W., Ruangsuj, P., Mangmool, S., Promptmas, C., & Yasawong, M. (2019). A Novel Nucleic Lateral Flow Assay for Screening phaR-Containing Bacillus spp. *Journal of Microbiology and Biotechnology*. <https://doi.org/10.4014/jmb.1907.07045>
4. Pecchia, S., & Da Lio, D. (2018). Development of a rapid PCR-Nucleic Acid Lateral Flow Immunoassay (PCR-NALFIA) based on rDNA IGS sequence analysis for the detection of *Macrophomina phaseolina* in soil. *Journal of Microbiological Methods*, 151, 118–128. <https://doi.org/10.1016/j.mimet.2018.06.010>
5. Zaky, W. I., Tomaino, F. R., Pilotte, N., Laney, S. J., & Williams, S. A. (2018). Backpack PCR: A point-of-collection diagnostic platform for the rapid detection of *Brugia* parasites in mosquitoes. *PLoS Neglected Tropical Diseases*, 12(11). <https://doi.org/10.1371/journal.pntd.0006962>
6. Janz, V., Schoon, J., Morgenstern, C., Preininger, B., Reinke, S., Duda, G., Breitbach, A., Perka, C. F., & Geissler, S. (2018). Rapid detection of periprosthetic joint infection using a combination of 16s rDNA polymerase chain reaction and lateral flow immunoassay: A Pilot Study. *Bone & joint research*, 7(1), 12–19. <https://doi.org/10.1302/2046-3758.71.BJR-2017-0103.R2>
7. Muangsuwan, W., Ruangsuj, P., Chaichanachaicharn, P., & Yasawong, M. (2015). A novel nucleic lateral flow assay for screening of PHA-producing haloarchaea. *Journal of Microbiological Methods*, 116, 8–14. <https://doi.org/10.1016/j.mimet.2015.06.012>
8. Wong, Y. L., Wong, K. L., & Shaw, P. C. (2015). Rapid authentication of *Cordyceps* by lateral flow dipstick. *Journal of Pharmaceutical and Biomedical Analysis*, 111, 306–310. <https://doi.org/10.1016/j.jpba.2015.04.003> **multiplex**
9. Breitbach A., Jacob F., Hutzler M., Koob J. (2017). Simple Rapid Test for Establishing the Presence of Beer-Spoilage Organisms. Brauwelt International. [Link to article](#). **multiplex**

Loop mediated isothermal amplification (LAMP)

1. Yongkiettrakul, S., Kolié, F. R., Kongkasuriyachai, D., Sattabongkot, J., Nguitragool, W., Nawattanapaibool, N., ... Buates, S. (2020). Validation of PfSNP-LAMP-Lateral Flow Dipstick for Detection of Single Nucleotide Polymorphism Associated with Pyrimethamine Resistance in *Plasmodium falciparum*. *Diagnostics*, 10(11), 948. <https://doi.org/10.3390/diagnostics10110948>
2. Ledlod, S., Bunroddith, K., Areekit, S., Santiwatanakul, S., & Chansiri, K. (2020). Development of a duplex lateral flow dipstick test for the detection and differentiation of *Listeria* spp. and *Listeria monocytogenes* in meat products based on loop-mediated isothermal amplification. *Journal of Chromatography B: Analytical Technologies in the Biomedical and Life Sciences*, 1139. <https://doi.org/10.1016/j.jchromb.2019.121834>
3. Huang, H. L., Gao, W. F., Zhu, P., Zhou, C. X., Qiao, L. L., Dang, C. Y., ... Yan, X. J. (2020). Molecular method for rapid detection of the red tide dinoflagellate *Karenia mikimotoi* in the coastal region of Xiangshan Bay, China. *Journal of Microbiological Methods*, 168. <https://doi.org/10.1016/j.mimet.2019.105801>
4. Allgöwer, S. M., Hartmann, C. A., & Holzhauser, T. (2020). The development of highly specific and sensitive primers for the detection of potentially allergenic soybean (*Glycine max*) using loop-mediated isothermal amplification combined with lateral flow dipstick (LAMP-LFD). *Foods*, 9(4). <https://doi.org/10.3390/foods9040423>
5. Ağel, E., Sağcan, H., Ceyhan, İ., & Durmaz, R. (2020). Optimization of isothermal amplification method for *Mycobacterium tuberculosis* detection and visualization method for fieldwork. *Turkish Journal of Medical Sciences*, 50(4), 1069–1075. <https://doi.org/10.3906/sag-1910-6>
6. Thongphueak, D., Chansiri, K., Sriyapai, T., Areekit, S., Santiwatanakul, S., & Wangroongsarb, P. (2019). Development of the rapid test kit for the identification of *Campylobacter* spp. Based on Loop-mediated Isothermal Amplification (LAMP) in combination with a Lateral Flow Dipstick (LFD) and Gold Nano-DNA Probe (AuNPs). *Science and Technology Asia*, 24(1), 63–71. <https://doi.org/10.14456/scitechasia.2019.7>
7. Sagcan, H., & Turgut Kara, N. (2019). Detection of Potato ring rot Pathogen *Clavibacter michiganensis* subsp. *sepedonicus* by Loop-mediated isothermal amplification (LAMP) assay. *Scientific Reports*, 9(1). <https://doi.org/10.1038/s41598-019-56680-9>
8. Liu, L., Xu, Y., Zhong, W., Li, L., Li, W., & Xiao, Q. (2019). Comparison of three terminal detection methods based on loop mediated isothermal amplification (Lamp) assay for spring viremia of carp virus (svcv). *Turkish Journal of Fisheries and Aquatic Sciences*, 19(9), 805–816. https://doi.org/10.4194/1303-2712-v19_9_09
9. Mamba, T. S., Mbae, C. K., Kinyua, J., Mulinge, E., Mburugu, G. N., & Njiru, Z. K. (2018). Lateral Flow Loop-Mediated Isothermal Amplification Test with Stem Primers: Detection of *Cryptosporidium* Species in Kenyan Children Presenting with Diarrhea. *Journal of Tropical Medicine*, 2018. <https://doi.org/10.1155/2018/7659730>
10. Lalle, M., Possenti, A., Dubey, J. P., & Pozio, E. (2018). Loop-Mediated Isothermal Amplification-Lateral-Flow Dipstick (LAMP-LFD) to detect *Toxoplasma gondii* oocyst in ready-to-eat salad. *Food Microbiology*, 70, 137–142. <https://doi.org/10.1016/j.fm.2017.10.001>

Loop mediated isothermal amplification (LAMP)

11. Kanchanaphum, P. (2018). Time course of detection of human Male DNA from stained blood sample on various surfaces by loop mediated isothermal amplification and polymerase chain reaction. *BioMed Research International*, 2018. <https://doi.org/10.1155/2018/2981862>
12. Kumvongpin, R., Jearanaikoon, P., Wilailuckana, C., Sae-Ung, N., Prasongdee, P., Daduang, S., ... Daduang, J. (2017). Detection assay for HPV16 and HPV18 by loop-mediated isothermal amplification with lateral flow dipstick tests. *Molecular Medicine Reports*, 15(5), 3203–3209. <https://doi.org/10.3892/mmr.2017.6370>
13. Yang, Y., Li, Q., Wang, S., Chen, X., & Du, A. (2016). Rapid and sensitive detection of Babesia bovis and Babesia bigemina by loop-mediated isothermal amplification combined with a lateral flow dipstick. *Veterinary Parasitology*, 219, 71–76. <https://doi.org/10.1016/j.vetpar.2016.02.004>
14. Nawattanapaiboon, K., Prombun, P., Santanirand, P., Vongsakulyanon, A., Sriksirin, T., Sutapun, B., & Kiatpathomchai, W. (2016). Hemoculture and Direct Sputum Detection of mecA-Mediated Methicillin-Resistant Staphylococcus aureus by Loop-Mediated Isothermal Amplification in Combination With a Lateral-Flow Dipstick. *Journal of Clinical Laboratory Analysis*, 30(5), 760–767. <https://doi.org/10.1002/jcla.21935>
15. Yun, X., Hao, Y., Peng, Z., Hailong, H., Qianqian, H., & Jian, Z. (2015). Repaid and accurate detection of Khuskia oryzae by loopmediated isothermal amplification combined with a lateral flow dipstick. *Chemical Engineering Transactions*, 46, 877–882. <https://doi.org/10.3303/CET1546147>
16. Plaon, S., Longyant, S., Sithigorngul, P., & Chaivisuthangkura, P. (2015). Rapid and sensitive detection of Vibrio alginolyticus by loop-mediated isothermal amplification combined with a lateral flow dipstick targeted to the rpoX gene. *Journal of Aquatic Animal Health*, 27(3), 156–163. <https://doi.org/10.1080/08997659.2015.1037468>
17. Kusumawati, A., Tampubolon, I. D., Hendarta, N. Y., Salasia, S. I. O., Wanahari, T. A., Mappakaya, B. A., & Hartati, S. (2015). Use of reverse transcription loop-mediated isothermal amplification combined with lateral flow dipstick for an easy and rapid detection of Jembrana disease virus. *VirusDisease*, 26(3), 189–195. <https://doi.org/10.1007/s13337-015-0277-5> RNA
18. Yongkiettrakul, S., Jaroenram, W., Arunrut, N., Chareanchim, W., Pannengpetch, S., Suebsing, R., ... Kongkasuriyachai, D. (2014). Application of loop-mediated isothermal amplification assay combined with lateral flow dipstick for detection of Plasmodium falciparum and Plasmodium vivax. *Parasitology International*, 63(6), 777–784. <https://doi.org/10.1016/j.parint.2014.06.004>
19. Sun, Y. L., Yen, C. H., & Tu, C. F. (2014). Visual detection of canine parvovirus based on loop-mediated isothermal amplification combined with enzyme-linked immunosorbent assay and with lateral flow dipstick. *Journal of Veterinary Medical Science*, 76(4), 509–516. <https://doi.org/10.1292/jvms.13-0448>
20. Rigano, L. A., Malamud, F., Orce, I. G., Filippone, M. P., Marano, M. R., Do Amaral, A. M., ... Vojnov, A. A. (2014). Rapid and sensitive detection of Candidatus Liberibacter asiaticus by loop mediated isothermal amplification combined with a lateral flow dipstick. *BMC Microbiology*, 14(1). <https://doi.org/10.1186/1471-2180-14-86>

Loop mediated isothermal amplification (LAMP)

21. Chowdry, V. K., Luo, Y., Widén, F., Qiu, H. J., Shan, H., Belák, S., & Liu, L. (2014). Development of a loop-mediated isothermal amplification assay combined with a lateral flow dipstick for rapid and simple detection of classical swine fever virus in the field. *Journal of Virological Methods*, 197, 14–18. <https://doi.org/10.1016/j.jviromet.2013.11.013> **RNA**
22. Vaagt, F., Haase, I., & Fischer, M. (2013). Loop-mediated isothermal amplification (LAMP)-based method for rapid mushroom species identification. *Journal of Agricultural and Food Chemistry*, 61(8), 1833–1840. <https://doi.org/10.1021/jf304824b>
23. Khunthong, S., Jaroenram, W., Arunrut, N., Suebsing, R., Mungsantisuk, I., & Kiatpathomchai, W. (2013, March). Rapid and sensitive detection of shrimp yellow head virus by loop-mediated isothermal amplification combined with a lateral flow dipstick. *Journal of Virological Methods*, Vol. 188, pp. 51–56. <https://doi.org/10.1016/j.jviromet.2012.11.041>
24. Kaewphinit, T., Arunrut, N., Kiatpathomchai, W., Santiwatanakul, S., Jaratsing, P., & Chansiri, K. (2013). Detection of Mycobacterium tuberculosis by using loop-mediated isothermal amplification combined with a lateral flow dipstick in clinical samples. *BioMed Research International*, 2013. <https://doi.org/10.1155/2013/926230>
25. Tsai, S. M., Liu, H. J., Shien, J. H., Lee, L. H., Chang, P. C., & Wang, C. Y. (2012). Rapid and sensitive detection of infectious bursal disease virus by reverse transcription loop-mediated isothermal amplification combined with a lateral flow dipstick. *Journal of Virological Methods*, 181(1), 117–124. <https://doi.org/10.1016/j.jviromet.2011.09.002> **RNA**
26. Surasilp, T., Longyant, S., Rukpratanporn, S., Sridulyakul, P., Sithigorngul, P., & Chaivisuthangkura, P. (2011). Rapid and sensitive detection of *Vibrio vulnificus* by loop-mediated isothermal amplification combined with lateral flow dipstick targeted to rpoS gene. *Molecular and Cellular Probes*, 25(4), 158–163. <https://doi.org/10.1016/j.mcp.2011.04.001>
27. Prompamorn, P., Sithigorngul, P., Rukpratanporn, S., Longyant, S., Sridulyakul, P., & Chaivisuthangkura, P. (2011). The development of loop-mediated isothermal amplification combined with lateral flow dipstick for detection of *Vibrio parahaemolyticus*. *Letters in Applied Microbiology*, 52(4), 344–351. <https://doi.org/10.1111/j.1472-765X.2011.03007.x>
28. Njiru, Z. K. (2011). Rapid and sensitive detection of human African trypanosomiasis by loop-mediated isothermal amplification combined with a lateral-flow dipstick. *Diagnostic Microbiology and Infectious Disease*, 69(2), 205–209. <https://doi.org/10.1016/j.diagmicrobio.2010.08.026>
29. Arunrut, N., Seetang-Nun, Y., Phromjai, J., Panphut, W., & Kiatpathomchai, W. (2011). Rapid and sensitive detection of Laem-Singh virus by reverse transcription loop-mediated isothermal amplification combined with a lateral flow dipstick. *Journal of Virological Methods*, 177(1), 71–74. <https://doi.org/10.1016/j.jviromet.2011.06.020>
30. Soliman, H., & El-Matbouli, M. (2010). Loop mediated isothermal amplification combined with nucleic acid lateral flow strip for diagnosis of cyprinid herpes virus-3. *Molecular and Cellular Probes*, 24(1), 38–43. <https://doi.org/10.1016/j.mcp.2009.09.002>

Loop mediated isothermal amplification (LAMP)

31. Puthawibool, T., Senapin, S., Flegel, T. W., & Kiatpathomchai Wansika, W. (2010). Rapid and sensitive detection of *Macrobrachium rosenbergii* nodavirus in giant freshwater prawns by reverse transcription loop-mediated isothermal amplification combined with a lateral flow dipstick. *Molecular and Cellular Probes*, 24(5), 244–249. <https://doi.org/10.1016/j.mcp.2010.07.003> **RNA**
32. Njiru, Z. K., Ouma, J. O., Enyaru, J. C., & Dargantes, A. P. (2010). Loop-mediated Isothermal Amplification (LAMP) test for detection of *Trypanosoma evansi* strain B. *Experimental Parasitology*, 125(3), 196–201. <https://doi.org/10.1016/j.exppara.2010.01.017>
33. Nimitphak, T., Meemetta, W., Arunrut, N., Senapin, S., & Kiatpathomchai, W. (2010). Rapid and sensitive detection of *Penaeus monodon* nucleopolyhedrovirus (PemoNPV) by loop-mediated isothermal amplification combined with a lateral-flow dipstick. *Molecular and Cellular Probes*, Vol. 24, pp. 1–5. <https://doi.org/10.1016/j.mcp.2009.09.004>
34. Ding, W. C., Chen, J., Shi, Y. H., Lu, X. J., & Li, M. Y. (2010). Rapid and sensitive detection of infectious spleen and kidney necrosis virus by loop-mediated isothermal amplification combined with a lateral flow dipstick. *Archives of Virology*, 155(3), 385–389. <https://doi.org/10.1007/s00705-010-0593-4>
35. Puthawibool, T., Senapin, S., Kiatpathomchai, W., & Flegel, T. W. (2009). Detection of shrimp infectious myonecrosis virus by reverse transcription loop-mediated isothermal amplification combined with a lateral flow dipstick. *Journal of Virological Methods*, 156(1–2), 27–31. <https://doi.org/10.1016/j.jviromet.2008.10.018> **RNA**
36. Kikuchi, T., Aikawa, T., Oeda, Y., Karim, N., & Kanzaki, N. (2009). A rapid and precise diagnostic method for detecting the pinewood nematode *Bursaphelenchus xylophilus* by loop-mediated isothermal amplification. *Phytopathology*, 99(12), 1365–1369. <https://doi.org/10.1094/PHYTO-99-12-1365>
37. Jaroenram, W., Kiatpathomchai, W., & Flegel, T. W. (2009, April). Rapid and sensitive detection of white spot syndrome virus by loop-mediated isothermal amplification combined with a lateral flow dipstick. *Molecular and Cellular Probes*, Vol. 23, pp. 65–70. <https://doi.org/10.1016/j.mcp.2008.12.003>
38. Andrade, T. P. D., & Lightner, D. V. (2009). Development of a method for the detection of infectious myonecrosis virus by reverse-transcription loop-mediated isothermal amplification and nucleic acid lateral flow hybrid assay. *Journal of Fish Diseases*, 32(11), 911–924. <https://doi.org/10.1111/j.1365-2761.2009.01072.x> **RNA**
39. Nimitphak, T., Kiatpathomchai, W., & Flegel, T. W. (2008). Shrimp hepatopancreatic parvovirus detection by combining loop-mediated isothermal amplification with a lateral flow dipstick. *Journal of Virological Methods*, 154(1–2), 56–60. <https://doi.org/10.1016/j.jviromet.2008.09.003>
40. Kiatpathomchai, W., Jaroenram, W., Arunrut, N., Jitrapakdee, S., & Flegel, T. W. (2008). Shrimp Taura syndrome virus detection by reverse transcription loop-mediated isothermal amplification combined with a lateral flow dipstick. *Journal of Virological Methods*, 153(2), 214–217. <https://doi.org/10.1016/j.jviromet.2008.06.025> **RNA**



Recombinase Polymerase Amplification (RPA)

1. Wang, J., Li, R., Sun, X., Liu, L., Hao, X., Wang, J., & Yuan, W. (2020). Development and validation of the isothermal recombinase polymerase amplification assays for rapid detection of *Mycoplasma ovipneumoniae* in sheep. *BMC Veterinary Research*, 16(1). <https://doi.org/10.1186/s12917-020-02387-3>
2. Srisrattakarn, A., Tippayawat, P., Chanawong, A., Tavichakorntrakool, R., Daduang, J., Wonglakorn, L., ... Lulitanond, A. (2020). Direct detection of methicillin-resistant in *Staphylococcus* spp. in positive blood culture by isothermal recombinase polymerase amplification combined with lateral flow dipstick assay. *World Journal of Microbiology and Biotechnology*, 36(11). <https://doi.org/10.1007/s11274-020-02938-8>
3. Li, Z., Torres, J. E. P., Goossens, J., Stijlemans, B., Sterckx, Y. G. J., & Magez, S. (2020). Development of a recombinase polymerase amplification lateral flow assay for the detection of active *Trypanosoma evansi* infections. *PLoS Neglected Tropical Diseases*, 14(2). <https://doi.org/10.1371/journal.pntd.0008044>
4. Zhai, Y., Ma, P., Fu, X., Zhang, L., Cui, P., Li, H., ... Yang, X. (2020). A recombinase polymerase amplification combined with lateral flow dipstick for rapid and specific detection of African swine fever virus. *Journal of Virological Methods*, 285. <https://doi.org/10.1016/j.jviromet.2020.113885>
5. Jia, T., Yu, Y., & Wang, Y. (2020). A recombinase polymerase amplification-based lateral flow strip assay for rapid detection of genogroup II noroviruses in the field. *Archives of Virology*, 165(12), 2767–2776. <https://doi.org/10.1007/s00705-020-04798-x> **RNA**
6. Gao, X., Liu, X., Zhang, Y., Wei, Y., & Wang, Y. (2020). Rapid and visual detection of porcine deltacoronavirus by recombinase polymerase amplification combined with a lateral flow dipstick. *BMC Veterinary Research*, 16(1). <https://doi.org/10.1186/s12917-020-02341-3>
7. Azinheiro, S., Carvalho, J., Prado, M., & Garrido-Maestu, A. (2020). Application of recombinase polymerase amplification with lateral flow for a naked-eye detection of *Listeria monocytogenes* on food processing surfaces. *Foods*, 9(9). <https://doi.org/10.3390/foods9091249>
8. Qian, J., Boswell, S. A., Chidley, C., Lu, Z. xiang, Pettit, M. E., Gaudio, B. L., ... Springer, M. (2020). An enhanced isothermal amplification assay for viral detection. *Nature Communications*, 11(1). <https://doi.org/10.1038/s41467-020-19258-y>
9. Dai, T., Wang, A., Yang, X., Yu, X., Tian, W., Xu, Y., & Hu, T. (2020). PHYCI_587572: An RxLR effector gene and new biomarker in a recombinase polymerase amplification assay for rapid detection of *Phytophthora cinnamomi*. *Forests*, 11(3). <https://doi.org/10.3390/f11030306>
10. Ravindran, V. B., Khallaf, B., Surapaneni, A., Crosbie, N. D., Soni, S. K., & Ball, A. S. (2020). Detection of Helminth Ova in Wastewater Using Recombinase Polymerase Amplification Coupled to Lateral Flow Strips. *Water*, 12(3), 691. <https://doi.org/10.3390/w12030691>
11. Saxena, A., Pal, V., Tripathi, N. K., & Goel, A. K. (2019). Development of a rapid and sensitive recombinase polymerase amplification-lateral flow assay for detection of *Burkholderia mallei*. *Transboundary and Emerging Diseases*, 66(2), 1016–1022. <https://doi.org/10.1111/tbed.13126>



Recombinase Polymerase Amplification (RPA)

12. Yu, J., Shen, D., Dai, T., Lu, X., Xu, H., & Dou, D. (2019). Rapid and equipment-free detection of *Phytophthora capsici* using lateral flow strip-based recombinase polymerase amplification assay. *Letters in Applied Microbiology*, 69(1), 64–70. <https://doi.org/10.1111/lam.13166>
13. Ma, Q., Yao, J., Yuan, S., Liu, H., Wei, N., Zhang, J., & Shan, W. (2019). Development of a lateral flow recombinase polymerase amplification assay for rapid and visual detection of *Cryptococcus neoformans*/*C. gattii* in cerebral spinal fluid. *BMC Infectious Diseases*, 19(1), 1–9. <https://doi.org/10.1186/s12879-019-3744-6>
14. Hu, S., Zhong, H., Huang, W., Zhan, W., Yang, X., Tang, B., ... Luo, M. (2019). Rapid and visual detection of Group B streptococcus using recombinase polymerase amplification combined with lateral flow strips. *Diagnostic Microbiology and Infectious Disease*, 93(1), 9–13. <https://doi.org/10.1016/j.diagmicrobio.2018.07.011>
15. Li, T. T., Wang, J. L., Zhang, N. Z., Li, W. H., Yan, H. Bin, Li, L., ... Fu, B. Q. (2019). Rapid and visual detection of *Trichinella* Spp. Using a lateral flow strip-based recombinase polymerase amplification (LF-RPA) assay. *Frontiers in Cellular and Infection Microbiology*, 9(JAN), 1–8. <https://doi.org/10.3389/fcimb.2019.00001>
16. Dai, T., Hu, T., Yang, X., Shen, D., Jiao, B., Tian, W., & Xu, Y. (2019). A recombinase polymerase amplification lateral flow dipstick assay for rapid detection of the quarantine citrus pathogen in China, *Phytophthora hibernalis*. *PeerJ*, 2019(11). <https://doi.org/10.7717/peerj.8083>
17. Dai, T., Yang, X., Hu, T., Jiao, B., Xu, Y., Zheng, X., & Shen, D. (2019). Comparative evaluation of a novel recombinase polymerase amplification-lateral flow dipstick (Rpa-lfd) assay, lamp, conventional pcr, and leaf-disc baiting methods for detection of *phytophthora sojae*. *Frontiers in Microbiology*, 10(AUG). <https://doi.org/10.3389/fmicb.2019.01884>
18. Zhao, G., He, H., & Wang, H. (2019). Use of a recombinase polymerase amplification commercial kit for rapid visual detection of *Pasteurella multocida*. *BMC Veterinary Research*, 15(1), 1–8. <https://doi.org/10.1186/s12917-019-1889-6>
19. Peng, Y., Zheng, X., Kan, B., Li, W., Zhang, W., Jiang, T., ... Qin, A. (2019). Rapid detection of *Burkholderia pseudomallei* with a lateral flow recombinase polymerase amplification assay. *PLoS ONE*, 14(7). <https://doi.org/10.1371/journal.pone.0213416>
20. Rames, E. K., & Macdonald, J. (2019). Rapid assessment of viral water quality using a novel recombinase polymerase amplification test for human adenovirus. *Applied Microbiology and Biotechnology*, 103(19), 8115–8125. <https://doi.org/10.1007/s00253-019-10077-w>
21. Chia, J. Y., Tay, Y. P., Poh, D., Tay, B. H., Koh, E., Teo, J., ... Lee, C. L. K. (2019, August 10). Development and evaluation of a rapid on-site water pathogen detection system for water quality monitoring. *BioRxiv*. <https://doi.org/10.1101/731778>
22. Liu, M. zhi, Han, X. hu, Yao, L. quan, Zhang, W. kui, Liu, B. shan, & Chen, Z. liang. (2019). Development and application of a simple recombinase polymerase amplification assay for rapid point-of-care detection of feline herpesvirus type 1. *Archives of Virology*, 164(1), 195–200. <https://doi.org/10.1007/s00705-018-4064-7>



Recombinase Polymerase Amplification (RPA)

23. Wu, L., Ye, L., Wang, Z., Cui, Y., & Wang, J. (2019). Utilization of recombinase polymerase amplification combined with a lateral flow strip for detection of *Perkinsus beihaiensis* in the oyster *Crassostrea hongkongensis*. *Parasites & Vectors*, 12(1), 360. <https://doi.org/10.1186/s13071-019-3624-3>
24. Zhang, T. T., Liu, M. Z., Yin, R. H., Yao, L. Q., Liu, B. S., & Chen, Z. L. (2019). Rapid and simple detection of *Glaesserella parasuis* in synovial fluid by recombinase polymerase amplification and lateral flow strip. *BMC Veterinary Research*, 15(1), 1–7. <https://doi.org/10.1186/s12917-019-2039-x>
25. Zhao, C., Sun, F., Li, X., Lan, Y., Du, L., Zhou, T., & Zhou, Y. (2019). Reverse transcription-recombinase polymerase amplification combined with lateral flow strip for detection of rice black-streaked dwarf virus in plants. *Journal of Virological Methods*, 263, 96–100. <https://doi.org/10.1016/j.jviromet.2018.11.001>
26. Zhao, G., Wang, H., Hou, P., He, C., & He, H. (2018). Rapid visual detection of *Mycobacterium avium* subsp. *Paratuberculosis* by recombinase polymerase amplification combined with a lateral flow dipstick. *Journal of Veterinary Science*, 19(2), 242–250. <https://doi.org/10.4142/jvs.2018.19.2.242>
27. Friedrich, R., Rappold, E., Bogdan, C., & Held, J. (2018). Comparative analysis of the wako B-glucan test and the fungitell assay for diagnosis of candidemia and *Pneumocystis jirovecii* pneumonia. *Journal of Clinical Microbiology*, 56(9). <https://doi.org/10.1128/JCM>
28. Gao, W., Huang, H., Zhu, P., Yan, X., Fan, J., Jiang, J., & Xu, J. (2018). Recombinase polymerase amplification combined with lateral flow dipstick for equipment-free detection of *Salmonella* in shellfish. *Bioprocess and Biosystems Engineering*, 41(5), 603–611. <https://doi.org/10.1007/s00449-018-1895-2>
29. Ahmed, F. A., Larrea-Sarmiento, A., Alvarez, A. M., & Arif, M. (2018). Genome-informed diagnostics for specific and rapid detection of *Pectobacterium* species using recombinase polymerase amplification coupled with a lateral flow device. *Scientific Reports*, 8(1). <https://doi.org/10.1038/s41598-018-34275-0>
30. Liu, Y., Wang, X. Y., Wei, X. M., Gao, Z. T., & Han, J. P. (2018). Rapid Authentication of *Ginkgo biloba* Herbal Products Using the Recombinase Polymerase Amplification Assay. *Scientific Reports*, 8(1). <https://doi.org/10.1038/s41598-018-26402-8>
31. Cui, J., Zhao, Y., Sun, Y., Yu, L., Liu, Q., Zhan, X., ... Zhao, J. (2018). Detection of *Babesia gibsoni* in dogs by combining recombinase polymerase amplification (RPA) with lateral flow (LF) dipstick. *Parasitology Research*, 117(12), 3945–3951. <https://doi.org/10.1007/s00436-018-6104-3>
32. Soliman, H., & El-Matbouli, M. (2018). Rapid detection and differentiation of carp oedema virus and cyprinid herpes virus-3 in koi and common carp. *Journal of Fish Diseases*, 41(5), 761–772. <https://doi.org/10.1111/jfd.12774> **multitplex**
33. Qi, Y., Yin, Q., Shao, Y., Cao, M., Li, S., Chen, H., ... Li, Y. (2018). Development of a rapid and visual nucleotide detection method for a Chinese epidemic strain of *Orientia tsutsugamushi* based on recombinase polymerase amplification assay and lateral flow test. *International Journal of Infectious Diseases*, 70, 42–50. <https://doi.org/10.1016/j.ijid.2018.03.003>

Recombinase Polymerase Amplification (RPA)

34. Zhao, G., Hou, P., Huan, Y., He, C., Wang, H., & He, H. (2018). Development of a recombinase polymerase amplification combined with a lateral flow dipstick assay for rapid detection of the *Mycoplasma bovis*. *BMC Veterinary Research*, 14(1). <https://doi.org/10.1186/s12917-018-1703-x> **RNA**
35. Lai, M. Y., Ooi, C. H., & Lau, Y. L. (2018). Recombinase polymerase amplification combined with a lateral flow strip for the detection of plasmodium knowlesi. *American Journal of Tropical Medicine and Hygiene*, 98(3), 700–703. <https://doi.org/10.4269/ajtmh.17-0738> **multiplex**
36. Wang, H., Hou, P., Zhao, G., Yu, L., Gao, Y. W., & He, H. (2018). Development and evaluation of serotype-specific recombinase polymerase amplification combined with lateral flow dipstick assays for the diagnosis of foot-and-mouth disease virus serotype A, O and Asia1. *BMC veterinary research*, 14(1), 359. <https://doi.org/10.1186/s12917-018-1644-4>
37. Hou, P., Zhao, G., Wang, H., He, C., & He, H. (2018). Rapid detection of bovine viral diarrhoea virus using recombinase polymerase amplification combined with lateral flow dipstick assays in bulk milk. *Vet. Arhiv*, 88(5), 627–642. <https://doi.org/10.24099/vet.arhiv.0145>
38. James, A. S., Todd, S., Pollak, N. M., Marsh, G. A., & Macdonald, J. (2018). Ebolavirus diagnosis made simple, comparable and faster than molecular detection methods: Preparing for the future. *Virology Journal*, 15(1). <https://doi.org/10.1186/s12985-018-0985-8> **RNA**
39. Qi, Y., Shao, Y., Rao, J., Shen, W., Yin, Q., Li, X., ... Li, Y. (2018). Development of a rapid and visual detection method for *Rickettsia rickettsii* combining recombinase polymerase assay with lateral flow test. *PLoS ONE*, 13(11). <https://doi.org/10.1371/journal.pone.0207811>
40. Soliman, H., Kumar, G., & El-Matbouli, M. (2018). Recombinase polymerase amplification assay combined with a lateral flow dipstick for rapid detection of *Tetracapsuloides bryosalmonae*, the causative agent of proliferative kidney disease in salmonids. *Parasites and Vectors*, 11(1), 1–8. <https://doi.org/10.1186/s13071-018-2825-5>
41. Wang, H., Sun, M., Xu, D., Podok, P., Xie, J., Jiang, Y., & Lu, L. (2018). Rapid visual detection of cyprinid herpesvirus 2 by recombinase polymerase amplification combined with a lateral flow dipstick. *Journal of Fish Diseases*, 41(8), 1201–1206. <https://doi.org/10.1111/jfd.12808>
42. Tu, P.-A. (2018). A Recombinase Polymerase Amplification Lateral Flow Dipstick for Field Diagnosis of Bovine Leukemia Virus Infection and its Effectiveness Compared to iiPCR. and ELISA. *Journal of Antivirals & Antiretrovirals*, 10(3). <https://doi.org/10.4172/1948-5964.1000178>
43. Qi, Y., Yin, Q., Shao, Y., Li, S., Chen, H., Shen, W., ... Li, Y. (2018). Rapid and Visual Detection of *Coxiella burnetii* Using Recombinase Polymerase Amplification Combined with Lateral Flow Strips. *BioMed Research International*, 2018. <https://doi.org/10.1155/2018/6417354>
44. Yang, Y., Qin, X., Song, Y., Zhang, W., Hu, G., Dou, Y., ... Zhang, Z. (2017). Development of real-time and lateral flow strip reverse transcription recombinase polymerase Amplification assays for rapid detection of peste des petits ruminants virus. *Virology Journal*, 14(1), 1–10. <https://doi.org/10.1186/s12985-017-0688-6> **RNA**



Recombinase Polymerase Amplification (RPA)

45. Hou, P., Wang, H., Zhao, G., He, C., & He, H. (2017). Rapid detection of infectious bovine Rhinotracheitis virus using recombinase polymerase amplification assays. *BMC Veterinary Research*, 13(1), 1–9. <https://doi.org/10.1186/s12917-017-1284-0>
46. Yang, Y., Qin, X., Sun, Y., Cong, G., Li, Y., & Zhang, Z. (2017). Development of Isothermal Recombinase Polymerase Amplification Assay for Rapid Detection of Porcine Circovirus Type 2. *BioMed Research International*, 2017. <https://doi.org/10.1155/2017/8403642>
47. Yang, Y., Qin, X., Zhang, X., Zhao, Z., Zhang, W., Zhu, X., ... Zhang, Z. (2017). Development of real-time and lateral flow dipstick recombinase polymerase amplification assays for rapid detection of goatpox virus and sheeppox virus. *Virology Journal*, 14(1), 1–8. <https://doi.org/10.1186/s12985-017-0792-7>
48. Szántó-Egész, R., Jánosi, A., Mohr, A., Szalai, G., Szabó, E. K., Micsinai, A., ... Zsolnai, A. (2016). Breed-Specific Detection of Mangalica Meat in Food Products. *Food Analytical Methods*, 9(4), 889–894. <https://doi.org/10.1007/s12161-015-0261-0>
49. Saldarriaga, O. A., Castellanos-Gonzalez, A., Porrozzi, R., Baldeviano, G. C., Lescano, A. G., de Los Santos, M. B., ... Travi, B. L. (2016). An Innovative Field-Applicable Molecular Test to Diagnose Cutaneous Leishmania Viannia spp. Infections. *PLoS Neglected Tropical Diseases*, 10(4). <https://doi.org/10.1371/journal.pntd.0004638>
50. Liu, W., Liu, H. X., Zhang, L., Hou, X. X., Wan, K. L., & Hao, Q. (2016). A novel isothermal assay of *Borrelia burgdorferi* by recombinase polymerase amplification with lateral flow detection. *International Journal of Molecular Sciences*, 17(8). <https://doi.org/10.3390/ijms17081250>
51. Yang, Y., Qin, X., Zhang, W., Li, Y., & Zhang, Z. (2016). Rapid and specific detection of porcine parvovirus by isothermal recombinase polymerase amplification assays. *Molecular and Cellular Probes*, 30(5), 300–305. <https://doi.org/10.1016/j.mcp.2016.08.011>
52. Yang, Y., Qin, X., Wang, G., Jin, J., Shang, Y., & Zhang, Z. (2016). Development of an isothermal amplification-based assay for rapid visual detection of an Orf virus. *Virology Journal*, 13(1), 1–6. <https://doi.org/10.1186/s12985-016-0502-x>
53. Wu, Y. D., Zhou, D. H., Zhang, L. X., Zheng, W. Bin, Ma, J. G., Wang, M., ... Xu, M. J. (2016). Recombinase polymerase amplification (RPA) combined with lateral flow (LF) strip for equipment-free detection of *Cryptosporidium* spp. oocysts in dairy cattle feces. *Parasitology Research*, 115(9), 3551–3555. <https://doi.org/10.1007/s00436-016-5120-4>
54. Sun, K., Xing, W., Yu, X., Fu, W., Wang, Y., Zou, M., ... Xu, D. (2016). Recombinase polymerase amplification combined with a lateral flow dipstick for rapid and visual detection of *Schistosoma japonicum*. *Parasites and Vectors*, 9(1), 1–9. <https://doi.org/10.1186/s13071-016-1745-5>
55. Rosser, A., Rollinson, D., Forrest, M., & Webster, B. L. (2015). Isothermal Recombinase Polymerase amplification (RPA) of *Schistosoma haematobium* DNA and oligochromatographic lateral flow detection. *Parasites and Vectors*, 8(1). <https://doi.org/10.1186/s13071-015-1055-3>



Recombinase Polymerase Amplification (RPA)

56. Chao, C. C., Belinskaya, T., Zhang, Z., & Ching, W. M. (2015). Development of recombinase polymerase amplification assays for detection of orientia tsutsugamushi or rickettsia typhi. *PLoS Neglected Tropical Diseases*, 9(7), 1–21. <https://doi.org/10.1371/journal.pntd.0003884>
57. Castellanos-Gonzalez, A., Saldarriaga, O. A., Tartaglino, L., Gacek, R., Temple, E., Sparks, H., ... Travi, B. L. (2015). A novel molecular test to diagnose canine visceral leishmaniasis at the point of care. *American Journal of Tropical Medicine and Hygiene*, 93(5), 970–975. <https://doi.org/10.4269/ajtmh.15-0145>
58. Cordray, M. S., & Richards-Kortum, R. R. (2015). A paper and plastic device for the combined isothermal amplification and lateral flow detection of Plasmodium DNA. *Malaria Journal*, 14(1). <https://doi.org/10.1186/s12936-015-0995-6>
59. Nair, G., Rebolledo, M., Clinton White, A., Crannell, Z., Rebecca Richards-Kortum, R., Elizabeth Pinilla, A., ... Castellanos-Gonzalez, A. (2015). Detection of entamoeba histolytica by recombinase polymerase amplification. *American Journal of Tropical Medicine and Hygiene*, 93(3), 591–595. <https://doi.org/10.4269/ajtmh.15-0276>
60. Crannell, Z. A., Rohrman, B., & Richards-Kortum, R. (2014). Equipment-free incubation of recombinase polymerase amplification reactions using body heat. *PLoS ONE*, 9(11). <https://doi.org/10.1371/journal.pone.0112146>
61. Crannell, Z. A., Castellanos-Gonzalez, A., Irani, A., Rohrman, B., White, A. C., & Richards-Kortum, R. (2014). Nucleic acid test to diagnose cryptosporidiosis: Lab assessment in animal and patient specimens. *Analytical Chemistry*, 86(5), 2565–2571. <https://doi.org/10.1021/ac403750z>
62. Escadafal, C., Faye, O., Sall, A. A., Faye, O., Weidmann, M., Strohmeier, O., ... Patel, P. (2014). Rapid Molecular Assays for the Detection of Yellow Fever Virus in Low-Resource Settings. *PLoS Neglected Tropical Diseases*, 8(3). <https://doi.org/10.1371/journal.pntd.0002730>
63. Kersting, S., Rausch, V., Bier, F. F., & Von Nickisch-Rosenegk, M. (2014). Rapid detection of Plasmodium falciparum with isothermal recombinase polymerase amplification and lateral flow analysis. *Malaria Journal*, 13(1). <https://doi.org/10.1186/1475-2875-13-99>
64. Rohrman, B. A., & Richards-Kortum, R. R. (2012). A paper and plastic device for performing recombinase polymerase amplification of HIV DNA. *Lab on a Chip*, 12(17), 3082–3088. <https://doi.org/10.1039/c2lc40423k>
65. Piepenburg, O., Williams, C. H., Stemple, D. L., & Armes, N. A. (2006). DNA detection using recombination proteins. *PLoS Biology*, 4(7), 1115–1121. <https://doi.org/10.1371/journal.pbio.0040204>



CRISPR/Cas-based Nuclease Assays

1. Ali, Z., Aman, R., Mahas, A., Rao, G. S., Tehseen, M., Marsic, T., ... Mahfouz, M. M. (2020). iSCAN: An RT-LAMP-coupled CRISPR-Cas12 module for rapid, sensitive detection of SARS-CoV-2. *Virus Research*, 288. <https://doi.org/10.1016/j.virusres.2020.19812> **Cas12 RNA**
2. Azhar, M., Phutela, R., Kumar, M., Hussain Ansari, A., Rauthan, R., Gulati, S., ... Maiti, S. (2020). Rapid, accurate, nucleobase detection using FnCas9 Single sentence summary A method to identify nucleotide sequence or nucleobase identity using FnCas9 and its implementation in the rapid and accurate diagnosis of SARS-CoV-2. *MedRxiv*, 2020.09.13.20193581. <https://doi.org/10.1101/2020.09.13.20193581> **fnCas9 RNA**
3. Barnes, K. G., Lachenauer, A. E., Nitido, A., Siddiqui, S., Gross, R., Beitzel, B., ... Sabeti, P. C. (2020). Deployable CRISPR-Cas13a diagnostic tools to detect and report Ebola and Lassa virus cases in real-time. *Nature Communications*, 11(1). <https://doi.org/10.1038/s41467-020-17994-9> **Cas13 RNA**
4. Brogan, D. J., Chaverra-Rodriguez, D., Lin, C. P., Smidler, A. L., Yang, T., Alcantara, L. M., ... Akbari, O. S. (2020). A Sensitive, Rapid, and Portable CasRx-based Diagnostic Assay for SARS-CoV-2. *MedRxiv: The Preprint Server for Health Sciences*. <https://doi.org/10.1101/2020.10.14.20212795> **Cas13d RNA**
5. Broughton, J. P., Deng, X., Yu, G., Fasching, C. L., Servellita, V., Singh, J., ... Chiu, C. Y. (2020). CRISPR-Cas12-based detection of SARS-CoV-2. *Nature Biotechnology*, 38(7), 870–874. <https://doi.org/10.1038/s41587-020-0513-4> **Cas12 RNA**
6. Curti, L., Pereyra-Bonnet, F., & Gimenez, C. A. (2020, March 2). An ultrasensitive, rapid, and portable coronavirus SARS-CoV-2 sequence detection method based on CRISPR-Cas12. *BioRxiv*. <https://doi.org/10.1101/2020.02.29.971127> **Cas12 RNA**
7. Hao, L., Zhao, R., Ngambenjawong, C., Fleming, H., & Bhatia, S. (2020). CRISPR-Cas-amplified urine biomarkers for multiplexed and portable cancer diagnostics. *BioRxiv*. <https://doi.org/10.1101/2020.06.17.157180> **Cas12**
8. Joung, J., Ladha, A., Saito, M., Kim, N.-G., Woolley, A. E., Segel, M., ... Zhang, F. (2020). Detection of SARS-CoV-2 with SHERLOCK One-Pot Testing. *New England Journal of Medicine*, 383(15), 1492–1494. <https://doi.org/10.1056/nejmc2026172> **Cas13 RNA**
9. Kaminski, M. M., Alcantar, M. A., Lape, I. T., Greensmith, R., Huske, A. C., Valeri, J. A., ... Collins, J. J. (2020). A CRISPR-based assay for the detection of opportunistic infections post-transplantation and for the monitoring of transplant rejection. *Nature Biomedical Engineering*, 4(6), 601–609. <https://doi.org/10.1038/s41551-020-0546-5> **Cas13**



CRISPR/Cas-based Nuclease Assays

10. Kanitchinda, S., Srisala, J., Suebsing, R., Prachumwat, A., & Chaijarasphong, T. (2020). CRISPR-Cas fluorescent cleavage assay coupled with recombinase polymerase amplification for sensitive and specific detection of *Enterocytozoon hepatopenaei*. *Biotechnology Reports*, 27. <https://doi.org/10.1016/j.btre.2020.e00485> **Cas12**
11. Lee, R. A., De Puig, H., Nguyen, P. Q., Angenent-Mari, N. M., Donghia, N. M., McGee, J. P., ... Collins, J. J. (2020). Ultrasensitive CRISPR-based diagnostic for field-applicable detection of *Plasmodium* species in symptomatic and asymptomatic malaria. *Proceedings of the National Academy of Sciences of the United States of America*, 117(41), 25722–25731. <https://doi.org/10.1073/pnas.2010196117> **Cas13**
12. Lu, S., Li, F., Chen, Q., Wu, J., Duan, J., Lei, X., ... Yin, H. (2020). Rapid detection of African swine fever virus using Cas12a-based portable paper diagnostics. *Cell Discovery*, 6(1). <https://doi.org/10.1038/s41421-020-0151-5>
13. Metsky, H. C., Freije, C. A., Kosoko-Thoroddsen, T. S. F., Sabeti, P. C., & Myhrvold, C. (2020). CRISPR-based surveillance for COVID-19 using genomically-comprehensive machine learning design. *BioRxiv*, 2020.02.26.967026. <https://doi.org/10.1101/2020.02.26.967026> **Cas12 RNA**
14. Nguyen, L. T., Smith, B. M., & Jain, P. K. (2020). Author Correction: Enhancement of trans-cleavage activity of Cas12a with engineered crRNA enables amplified nucleic acid detection (Nature Communications, (2020), 11, 1, (4906), 10.1038/s41467-020-18615-1). *Nature Communications*, 11(1). <https://doi.org/10.1038/s41467-020-20117-z> **Cas12**
15. Ooi, K. H., Tay, J. W. D., Teo, S. Y., Liu, M. M., Kaewsapsak, P., Jin, S., ... Tan, M. H. (2020). A CRISPR-based SARS-CoV-2 diagnostic assay that is robust against viral evolution and RNA editing. *BioRxiv*. <https://doi.org/10.1101/2020.07.03.185850> **Cas12**
16. Patchesung, M., Jantarug, K., Pattama, A., Aphicho, K., Suraritdechachai, S., Meesawat, P., ... Uttamapinant, C. (2020). Clinical validation of a Cas13-based assay for the detection of SARS-CoV-2 RNA. *Nature Biomedical Engineering*, 4(12), 1140–1149. <https://doi.org/10.1038/s41551-020-00603-x> **Cas13**
17. Rauch, J. N., Valois, E., Solley, S. C., Braig, F., Lach, R. S., Audouard, M., ... Wilson, M. Z. (2020). A Scalable, easy-to-deploy, protocol for Cas13-based detection of SARS-CoV-2 genetic material. *BioRxiv*. <https://doi.org/10.1101/2020.04.20.052159> **Cas13 RNA**
18. Schermer, B., Fabretti, F., Damagnez, M., Cristanziano, V. Di, Heger, E., Arjune, S., ... Müller, R. U. (2020). Rapid SARS-CoV-2 testing in primary material based on a novel multiplex LAMP assay. *MedRxiv*. <https://doi.org/10.1101/2020.06.18.20130377> **RNA**
19. Wang, X., Ji, P., Fan, H., Dang, L., Wan, W., Liu, S., ... Liao, M. (2020). CRISPR/Cas12a technology combined with immunochromatographic strips for portable detection of African swine fever virus. *Communications Biology*, 3(1). <https://doi.org/10.1038/s42003-020-0796-5> **Cas12**



CRISPR/Cas-based Nuclease Assays

20. Wang, X., Xiong, E., Tian, T., Cheng, M., Lin, W., Wang, H., ... Zhou, X. (2020). Clustered Regularly Interspaced Short Palindromic Repeats/Cas9-Mediated Lateral Flow Nucleic Acid Assay. *ACS Nano*, 14(2), 2497–2508. <https://doi.org/10.1021/acsnano.0c00022> **Cas9**
21. Yoshimi, K., Takeshita, K., Yamayoshi, S., Shibumura, S., Yamauchi, Y., Yamamoto, M., ... Mashimo, T. (2020). Rapid and accurate detection of novel coronavirus SARS-CoV-2 using CRISPR-Cas3. *MedRxiv*. <https://doi.org/10.1101/2020.06.02.20119875> **Cas3**
22. Abudayyeh, O. O., Gootenberg, J. S., Kellner, M. J., & Zhang, F. (2019). Nucleic Acid Detection of Plant Genes Using CRISPR-Cas13. *The CRISPR Journal*, 2(3), 165–171. <https://doi.org/10.1089/crispr.2019.0011> **Cas13**
23. Bai, J., Lin, H., Li, H., Zhou, Y., Liu, J., Zhong, G., ... Huang, L. (2019). Cas12a-Based On-Site and Rapid Nucleic Acid Detection of African Swine Fever. *Frontiers in Microbiology*, 10. <https://doi.org/10.3389/fmicb.2019.02830> **Cas12**
24. Kellner, M. J., Koob, J. G., Gootenberg, J. S., Abudayyeh, O. O., & Zhang, F. (2019). SHERLOCK: nucleic acid detection with CRISPR nucleases. *Nature Protocols*, 14(10), 2986–3012. <https://doi.org/10.1038/s41596-019-0210-2> **Cas13**
25. Sullivan, T. J., Dhar, A. K., Cruz-Flores, R., & Bodnar, A. G. (2019). Rapid, CRISPR-Based, Field-Deployable Detection Of White Spot Syndrome Virus In Shrimp. In *Scientific Reports* (Vol. 9). <https://doi.org/10.1038/s41598-019-56170-y> **Cas13**
26. Tsou, J. H., Leng, Q., & Jiang, F. (2019). A CRISPR Test for Detection of Circulating Nuclei Acids. *Translational Oncology*, 12(12), 1566–1573. <https://doi.org/10.1016/j.tranon.2019.08.011> **Cas12**
27. Gootenberg, J. S., Abudayyeh, O. O., Kellner, M. J., Joung, J., Collins, J. J., & Zhang, F. (2018). Multiplexed and portable nucleic acid detection platform with Cas13, Cas12a and Csm6. *Science*, 360(6387), 439–444. <https://doi.org/10.1126/science.aag0179> **Cas13**
28. Myhrvold, C., Freije, C. A., Gootenberg, J. S., Abudayyeh, O. O., Metsky, H. C., Durbin, A. F., ... Sabeti, P. C. (2018). Field-deployable viral diagnostics using CRISPR-Cas13. In *Science* (Vol. 360). <https://doi.org/10.1126/science.aas8836> **Cas13**
29. Gootenberg, J. S., Abudayyeh, O. O., Lee, J. W., Essletzbichler, P., Dy, A. J., Joung, J., ... Zhang, F. (2017). Nucleic acid detection with CRISPR-Cas13a/C2c2. In *Science* (Vol. 356). <https://doi.org/10.1126/science.aam9321> **Cas13**

Protocols for the detection of SARS-CoV-2 with SHERLOCK or DETECTR

30. Zhang F., Abudayyeh O. O., Gootenberg J.S., Mathers L. (2020). A protocol for detection of COVID-19 using CRISPR diagnostics. Broad Communications. [Link to protocol](#). **Cas13**
31. Broughton J. P., Wayne D., Fasching C. L., Singh J., Charles Y., Chen J. S. (2020). A protocol for rapid detection of the 2019 novel coronavirus SARS-CoV-2 using CRISPR diagnostics : SARS-CoV-2 DETECTR. MammothBiosciences. [Link to protocol](#). **Cas12**



Special Applications

Rolling Circle Amplification and Lateral Flow

1. Zingg, J. M., & Daunert, S. (2018). Trinucleotide rolling circle amplification: A novel method for the detection of RNA and DNA. *Methods and Protocols*, 1(2), 1–20. <https://doi.org/10.3390/mps1020015>

Aptamers and Lateral Flow

2. Fischer, C., Wessels, H., Paschke-Kratzin, A., & Fischer, M. (2017). Aptamers: Universal capture units for lateral flow applications. *Analytical Biochemistry*, 522, 53–60. <https://doi.org/10.1016/j.ab.2017.01.012>

Thermophilic helicase dependent isothermal Amplification (tHDA)

3. Zasada, A. A., Zacharczuk, K., Formińska, K., Wiatrzyk, A., Ziółkowski, R., & Malinowska, E. (2018). Isothermal DNA amplification combined with lateral flow dipsticks for detection of biothreat agents. *Analytical biochemistry*, 560, 60–66. <https://doi.org/10.1016/j.ab.2018.09.008>